
ASEM MEDICAL



Akatalug n yibenk n wallen

63 n yibenk

Afran azayaz n yibenk ASEM. Tulmisin d tsuddisin yellan ilaq ad ttwasentmen ɣur umselɣu. Akatalug-a ur yettbeddil ara isemliyen n umselɣu.

Yettwarna ass n 2026-11-10

Huvitz | ASEM

OCTavius



Semmus n yikayaden deg yiwet n tɣawsa n tugniwin n wallen.

Talɣut yettwaseknen

Tugniwin s waṭas n tarrayin

OCT, retinografi, angiyografi, biyumetri d topografi.

huvitz.com · 2026-10-05

Anṭaf OCT arurad

80 000 A-scans deg tasint i unṭaf OCT.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HBM-1



Biyumetri taḍwant d topografi n lqerniya i ukayad send agradd.

Talɣut yettwaseknen

Biyumetri taḍwant

Akat n teɣzi ɣef ugellus seg 14 ar 40 mm.

www.huvitz.com · 2026-10-05

Topografi Placido

24 n tneɣrutin i tkarḍa n lqerniya.

www.huvitz.com · 2026-10-05

Asiḍen n IOL

Asiḍen n IOL d tesleḍt Zernike deg yiwen.

www.huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HRK-9000A



Rifraksyun d ukayad n tmaẓirt n tiṭ deg yiwen n yibenk.

Talɣut yettwaseknen

Rifraksyun tawurmant

SPH seg −30 ar +25 D; CYL alamma 12 D.

huvitz.com · 2026-10-05

Tamaẓirt n tiṭ

Meybugrafi d tesleḍt n wakud n trẓi n tɣellayt n imeṭṭawen.

huvitz.com · 2026-10-05

Keratumetri n yiri

Akat n ukrif n lqerniya deg tlemmast d yiri.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Haag-Streit | ASEM

Eyestar 900



Biyumetri taḍwant d keratumetri s unagraw EyeSuite.

Talɣut yettwaseknen

Biyumetri OCT swept-source

Biyumetri OCT swept-source seg lqerniya ar retina.

haag-streit.com · 2026-10-17

Takarḍa n lqerniya

Allalen n topografi d pakimetri i usenqed n lqerniya.

haag-streit.com · 2026-10-17

Anṭaf awurman n yisefka

Anṭaf awurman yettusebda i ukayad n snat n wallen.

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Haag-Streit | ASEM

Elara 900



Taftilt s tecraḍ, s umikruskuḍ agalili d 5 n yisemeɣriyen.

Talɣut yettwaseknen

Tafat s usnay n tecraḍ

Tafat s usnay i tmugli n tɣessiwin n tiṭ.

us.haag-streit.com · 2026-10-17

Asenqed s uḍaw infraruz

Aḍaw infraruz i usenqed n tɣuddanin Meibomius.

us.haag-streit.com · 2026-10-17

Asenqed s umutur

Asenqed s umutur n teɣzi n yibenk d usemɣer.

us.haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · us.haag-streit.com

F.I.S.O. | ASEM

Easyvisus RM64 — 3 n yibenk



Tayunt n rifraksyun i 3 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 3 n yibenk

Tafelwit tettusebda i 3 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Asebded elektrik

Aseggem elektrik n tfelwit d ukersi.

www.fiso.info · 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

GMV | ASEM

Plexr® Plus



Tetiknuluji n plasma s 3 n yifassen i yiswiren yemgaraden.

Talɣut yettwaseknen

Tetiknuluji n plasma

Tazmert tettunefk s uqaws n plasma yettulemmesen.

gruppogmv.com · 2026-10-17

3 n yifassen

Kraḍ n yifassen: amellal, azegzaw d uzeggaɣ.

gruppogmv.com · 2026-10-17

Aseffeɣ war anermis

Aseffeɣ n tzemmert war anermas srid ɣef uglim.

gruppogmv.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · gruppogmv.com

BVI Medical | ASEM

R-Evo Smart® CR



Fakumulsifikasyun d fitriktumi n zdat d deffir.

Talɣut yettwaseknen

Agzay yedduklen

Fako d fitriktumi n zdat d deffir.

www.bvimedical.com · 2026-10-17

Sin n yizamulen n ussuḍeḍ

Tafrant gar tpumpi peristaltik d Venturi.

www.bvimedical.com · 2026-10-17

Fako s ultrasons

Aseffeɣ ameɣlal, s yinubag neɣ s ikellafen ɣer 40 kHz.

www.bvimedical.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.bvimedical.com

Haag-Streit | ASEM

METIS 900



Amikruskuḍ n ugzay n wallen s zoom d ulemmes s umutur.

Talɣut yettwaseknen

Reflex azeggaɣ coaxial

Tafat n reflex azeggaɣ coaxial i tmugli n ugzay.

us.haag-streit.com · 2026-10-17

Aḍran n tmugli s umutur

Aḍran s umutur i useggem n teɣmert n tmugli n ugzay.

us.haag-streit.com · 2026-10-17

Asenqed s ufus d uḍar

Asenqed n zoom d ulemmes s ufus neɣ s uḍar.

Ilmend n twila

us.haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · us.haag-streit.com

F.I.S.O. | ASEM

Easyvisus MS33 — 3 n yibenk



Tayunt n rifraksyun i 3 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 3 n yibenk

Tafelwit tettusebda i 3 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Asebded elektrik

Aseggem elektrik n tfelwit d ukersi.

www.fiso.info · 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Easyvisus Curve



Tayunt n rifraksyun i 2 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Tifranin n tfelwit

Lqimat yezziyen neɣ yettuzzuɣuren i 2 n yibenk.

Ilmend n twila

www.fiso.info · 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

GMV | ASEM

Plexr® PRO



Tetiknuluji n plasma s 2 n yifassen i yiswiren yemgaraden.

Talɣut yettwaseknen

Tetiknuluji n plasma

Tazmert tettunefk s uqaws n plasma yettulemmesen.

gruppogmv.com · 2026-10-17

2 n yifassen

Sin n yifassen: P One d P Two.

gruppogmv.com · 2026-10-17

Aseffeɣ war anermis

Aseffeɣ n tzemmert war anermas srid ɣef uglim.

gruppogmv.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · gruppogmv.com

PIROP A+B+CCT



Biyumetri, ultrasons d pakimetri deg yiwen n yibenk n wallen.

Talɣut yettwaseknen

Biyumetri s uskar A

Teɣzi ɣef ugellus d usiḍen n tzemmert n lentiy IOL.

echoson.pl · 2026-10-17

Tugniwin ultrasons s uskar B

Tugniwin ultrasons n tɣessiwin yellan daxel n tiṭ.

echoson.pl · 2026-10-17

Pakimetri CCT

Akat s ultrasons n tḍuḍi n lqerniya.

echoson.pl · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · echoson.pl

Huvitz | ASEM

HTR-1A



Rifraksyun, lqerniya d tḍeɣṭ deg tiṭ deg yiwet n tɣawsa.

Talɣut yettwaseknen

Rifraksyun d keratumetri

Sin n yikaten n ukayad n rifraksyun mlalen.

huvitz.com · 2026-10-05

Tonumetri war anermis

Akat n tḍeɣṭ deg tiṭ war anermis.

huvitz.com · 2026-10-05

Pakimetri n lqerniya

Akat n tḍuḍi n lqerniya deg yiwen n yibenk.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

HOCT-1 / HOCT-1F



OCT aspektral d retinografi war asemɣer n tberract deg yiwen.

Talɣut yettwaseknen

OCT aspektral

Tafrayt ɣef ugellus n 7 µm i tugniwin OCT.

huvitz.com · 2026-10-05

Retinografi deg yiwen

Tugna n retina war asemɣer n tberract.

HOCT-1F kan

huvitz.com · 2026-10-05

Asefrak n tugniwin

Asefrak n yikayaden s useɣẓan HIIS-1.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HFC-1



Retinografi war asemɣer n tberract s unṭaf apanurami.

Talɣut yettwaseknen

Tugniwin war asemɣer

Anṭaf n tugna n retina s wurti n 45°.

www.huvitz.com · 2026-10-05

Tamugli tapanuramit

Asdukel apanurami n tugniwin n retina.

www.huvitz.com · 2026-10-05

Taberract tameẓẓyant

Anṭaf seg ukatar n tberract n 3,3 mm.

www.huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HDR-9000



Arifraktur umḍin i ukayad asaguran d usiweḍ n yikaten.

Talɣut yettwaseknen

Rifraksyun tasagurant

SPH seg −29 ar +26,75 D i ukayad n tmugli.

huvitz.com · 2026-10-05

Aseɣti aprizmatik

Tazmert taprizmatikt alamma 20 Δ.

huvitz.com · 2026-10-05

Tuqqna war tineɣrist

Asemmeskel n yisefka s Wi-Fi neɣ RS-232.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HVS-1



Asenqed n rifraksyun n snat n wallen seg lbeɛd.

Talɣut yettwaseknen

Asenqed n snat n wallen

Akat n rifraksyun n snat n wallen seg lbeɛd.

huvitz.com · 2026-10-05

Azilal n rifraksyun

SPH seg −8 ar +8 D; CYL alamma 3 D.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HS-7500



Taḍwant tagalilit d semmus n yisemeɣriyen i ukayad s tecraḍ.

Talɣut yettwaseknen

Semmus n yisemeɣriyen

6×, 10×, 16×, 25×, 40× i usmṣada n tmugli.

huvitz.com · 2026-10-04

Tacraḍ n ukayad

Aldi n tecraḍ alamma 14 × 14 mm.

huvitz.com · 2026-10-04

Tistayin n tmugli

Tastayt tazegzawt cobalt d tazegzawt i ukayad.

huvitz.com · 2026-10-04

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

HRK-8100A



Rifraksyun tawurmant s keratumetri n tlemmast d yiri.

Talɣut yettwaseknen

Rifraksyun tawurmant

SPH seg −30 ar +25 D; CYL alamma 12 D.

huvitz.com · 2026-10-05

Lqerniya n yiri

Keratumetri n tlemmast d yiri deg yiwen n ukayad.

huvitz.com · 2026-10-05

Akrif n lentiyat

Akat n ukrif n lentiyat n unermis.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Haag-Streit | ASEM

Lenstar 900



Abiyumetr aḍwan s keratumetri d useɣẓan EyeSuite.

Talɣut yettwaseknen

Biyumetri taḍwant OLCR

Ikayaden ɣef ugellus n tiṭ i uheyyi n lentiy IOL.

haag-streit.com · 2026-10-17

Ikayaden n 2 n temnaḍin

Keratumetri n 2 n temnaḍin i usenqed n ukrif n Iqerniya.

haag-streit.com · 2026-10-17

Aheyyi IOL s EyeSuite

Asiḍen d tefrant n IOL deg EyeSuite IOL.

Ilmend n useggem n useɣẓan

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Haag-Streit | ASEM

Refractor 900



Arifraktur i ukayad n rifraksyun d tuqqna n yisefka.

Talɣut yettwaseknen

Rifraksyun tasagurant

Aseggem n lentiyat i yikayaden n rifraksyun.

haag-streit.com · 2026-10-17

Asenqed s ugdil anermas

Asenqed n ukayad s tfelwit n ugdil anermas.

haag-streit.com · 2026-10-17

Asemmeskel n yisefka

Asemmeskel n yisefka akked yibenkan d yikaramen n imuḍan.

Senqed amṣada n unagraw

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Haag-Streit | ASEM

Octopus 900



Perimetri ur nettbeddil ara; tiwuriwin-nniḍen ilmend n twila.

Talɣut yettwaseknen

Perimetri n wurti ummid

Ikayaden n wurti n tlemmast d yiri s yinagrawen usbiḍen.

haag-streit.com · 2026-10-17

Akayad n wurti s umussu

Ikayaden s umussu ilmend n lqem neɣ turagt.

Tafrant ilmend n lqem

haag-streit.com · 2026-10-17

Tasleḍt tamatut n unerni

Tasleḍt n unerni s yikayaden imseḍfaren n wurti n tmugli.

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Haag-Streit | ASEM

Octopus 600



Aperimetr s perimetri tusbiḍt SAP d perimetri Pulsar.

Talɣut yettwaseknen

Perimetri tusbiḍt n tlemmast

Ikayaden n wurti n tlemmast, amellal ɣef umellal.

haag-streit.com · 2026-10-17

Tasleḍt tamatut n unerni

Asenqed n yikayaden imseḍfaren d unerni-nsen amatu.

haag-streit.com · 2026-10-17

War tɣerɣert taberkant

Ikayaden zemren ad ilin berra n tɣerɣert taberkant.

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Haag-Streit | ASEM

BQ 900



Taftilt s tecraḍ, amikruskuḍ agalili, 5 n yisemeɣriyen d LED.

Talɣut yettwaseknen

Taḍwant n umikruskuḍ agalili

Semmus n yisemeɣriyen seg 6,3× ar 40×.

haag-streit.com · 2026-10-17

Tafat LED n tecraḍ

Tafat LED i tmugli s teftilt s tecraḍ.

Lqem LED

de.haag-streit.com · 2026-10-17

Tifranin n tugniwin

Anṭaf n tugniwin s uzegrir IM600 neɣ IM910.

Azegrir n tugniwin d tafrant

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.haag-streit.com

Haag-Streit | ASEM

BI 900



Taftilt s tecraḍ, amikruskuḍ Kepler, 2 n yisemeɣriyen d LED.

Talɣut yettwaseknen

Taḍwant n umikruskuḍ Kepler

Sin n yisemeɣriyen 10× d 16× i tmugli s tecraḍ.

us.haag-streit.com · 2026-10-17

Tastayt n umgired deg yiwen

Tastayt n umgired deg yiwen i yikayaden n yal ass.

us.haag-streit.com · 2026-10-17

Tifranin n tugniwin

Asedles n yikayaden s BI900 Imaging Set.

Tafeggagt n tugniwin d tafrant

us.haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.haag-streit.com

Haag-Streit | ASEM

Imaging Module 910



Azegrir n tugniwin i teftilt yemṣadan; twila ad ttwasenqed.

Talɣut yettwaseknen

Anṭaf n tugniwin s tallalt

Asenqed awurman n tafat d uldi i tallalt n unṭaf.

Taftilt yemṣadan ilaq

haag-streit.com · 2026-10-17

Asenqed n tkamyut iqerben

Isenqaden n tkamyut d tafat ɣer tama n joystick.

haag-streit.com · 2026-10-17

Sin n yiskaren n usdukel

Amahil iman-is neɣ asdukel daxel n EyeSuite.

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

F.I.S.O. | ASEM

Easyvisus RM64 — 2 n yibenk



Tayunt n rifraksyun i 2 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Asebded elektrik

Aseggem elektrik n tfelwit d ukersi.

www.fiso.info · 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Easyvisus RM64 Compact — 3 n yibenk



Tayunt n rifraksyun i 3 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 3 n yibenk

Tafelwit tettusebda i 3 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Anermis s ukersi s tirwiḍin

Tafelwit tettildi i ukayad n umuḍin ɣef ukersi s tirwiḍin.

www.fiso.info · 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

Easyvisus RM64 Compact — 2 n yibenk



Tayunt n rifraksyun i 2 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Anermis s ukersi s tirwiḍin

Tafelwit tettildi i ukayad n umuḍin ɣef ukersi s tirwiḍin.

www.fiso.info · 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Easyvisus MS33 — 2 n yibenk



Tayunt n rifraksyun i 2 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Asebded elektrik

Aseggem elektrik n tfelwit d ukersi.

www.fiso.info · 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Easyvisus One



Tayunt n rifraksyun i yiwen n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i yiwen n yibenk

Tafelwit tettusebda i yiwen n yibenk n ukayad.

www.fiso.info · 2026-10-17

Tafelwit s teɣzi ur nettbeddil ara

Aɣil s usewḥel amikanik deg yal amḍiq.

www.fiso.info · 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Easyvisus Two



Tayunt n rifraksyun i 2 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Alemmes n yibenk

Tafelwit s talɣa V d tarrayt n ulemmes n yibenk.

www.fiso.info · 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Panvisus Trio



Tayunt n rifraksyun i 3 n yibenk n ukayad.

Talɣut yettwaseknen

Amkan i 3 n yibenk

Tafelwit tettusebda i 3 n yibenk n ukayad.

www.fiso.info · 2026-10-17

Asebded elektrik

Aseggem elektrik n tfelwit d ukersi.

www.fiso.info · 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

www.fiso.info · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

BVI Medical | ASEM

R-Evo Smart® S



Fakumulsifikasyun d fitriktumi n zdat.

Talɣut yettwaseknen

Aferdis n zdat

Fakumulsifikasyun d fitriktumi n zdat.

www.bvimedical.com · 2026-10-17

Sin n yizamulen n ussuḍeḍ

Tafrant gar tpumpi peristaltik d Venturi.

www.bvimedical.com · 2026-10-17

Fako s ultrasons

Aseffeɣ ameɣlal, s yinubag neɣ s ikellafen ɣer 40 kHz.

www.bvimedical.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.bvimedical.com

BVI Medical | ASEM

R-Evo Smart® E



Fakumulsifikasyun d fitriktumi n zdat.

Talɣut yettwaseknen

Aferdis n zdat

Fakumulsifikasyun d fitriktumi n zdat.

www.bvimedical.com · 2026-10-17

Assuḍeḍ peristaltik

Assuḍeḍ s tpumpit peristaltik.

www.bvimedical.com · 2026-10-17

Fako s ultrasons

Aseffeɣ ameɣlal, s yinubag neɣ s ikellafen ɣer 40 kHz.

www.bvimedical.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.bvimedical.com

BVI Medical | ASEM

R-Evolution CR



Fakumulsifikasyun d fitriktumi n zdat d deffir.

Talɣut yettwaseknen

Agzay yedduklen

Fako d fitriktumi n zdat d deffir.

www.bvimedical.com · 2026-10-17

Sin n yizamulen n ussuḍeḍ

Tafrant gar tpumpi peristaltik d Venturi.

www.bvimedical.com · 2026-10-17

Fako s ultrasons

Aseffeɣ ameɣlal, s yinubag neɣ s ikellafen ɣer 40 kHz.

www.bvimedical.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.bvimedical.com

BVI Medical | ASEM

R-Evolution



Tanagrawt n fakumulsifikasyun s twuriwin n ugzay yettuseggamen.

Talɣut yettwaseknen

Asefrak n waman

Tafrant n ussuḍeḍ s uḍfar neɣ s tḍeɣṭ.

Ilmend n twila

www.bvimedical.com · 2026-10-17

Aseffeɣ n ultrasons

Iskaren s yinubag d yikellafen i fakumulsifikasyun.

Ilmend n twila

www.bvimedical.com · 2026-10-17

Tafat n ugzay

Kraḍ n tuffɣiwin LED ilelliyen s tistayin n yini.

Ilmend n twila

www.bvimedical.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · www.bvimedical.com

F.I.S.O. | ASEM

Table Hercule



Nermes ASEM i yisallen yettwasenqden ɣef temselɣut-a.

Tulmisin titiknikanin ilaq ad ttwasentmen.

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

Trial Lenses set with Metal Rim in Beech Wood Case



236 n lentiyat n ukayad n rifraksyun deg usenduq n usɣar.

Talɣut yettwaseknen

Rifraksyun tasagurant

Tafeggagt n lentiyat i usdukel n useɣti n tmugli.

www.fiso.info · 2026-10-05

236 n lentiyat n ukayad

Tafeggagt n 236 n lentiyat deg usenduq n usɣar.

www.fiso.info · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

51-B – 236 B Ensembles de lentilles d'essai



236 n lentiyat n ukayad n rifraksyun deg usenduq ABS.

Talɣut yettwaseknen

Rifraksyun tasagurant

Tafeggagt n lentiyat i usdukel n useɣti n tmugli.

www.fiso.info · 2026-10-05

236 n lentiyat n ukayad

Tafeggagt n 236 n lentiyat deg usenduq ABS.

www.fiso.info · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

lit chirurgical 77FY



Akersi n ugzay s yiferdisen yemmadawen i ugzay n wass.

Talɣut yettwaseknen

Asers i ugzay

Aseggem n yiri, n yiḍarren d uḍran n ukersi.

www.frastema.com · 2026-10-05

Amussu d usewḥel

Tirwiḍin yettunadin d usewḥel s uḍar.

www.frastema.com · 2026-10-05

Tifranin n usečči elektrik

Lqem s trisiti neɣ s tbatrit n 24 V, ilmend n twila.

Ilmend n twila

www.frastema.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

F.I.S.O. | ASEM

92-FA



Taburett s rrwaḍi d ukersi yettuseggmen i umkan n usenqed.

Talɣut yettwaseknen

Teɣzi tapnumatik

Amekanizm apnumatik iseggem teɣzi n ukersi seg 52 ar 72 cm.

Taburett FISO 92-FA. Aɣummu n ukersi d wayen ara d-yettunefken deg usiweḍ ilaq ad ttwasentmen.

fiso.info · 2026-10-18

Akersi imdewwer s uɣummu

Akersi imdewwer d ileqqaq, akat-is 33 cm.

fiso.info · 2026-10-18

Azadur s semmus n yifassen

Azadur s uɣummu krumi, akat-is 50 cm, s rrwaḍi.

fiso.info · 2026-10-18

Wali talqayt

Wali taɣbalut n umfaras · fiso.info

F.I.S.O. | ASEM

Panvisus Rondo



Tayunt Panvisus Rondo i sin n yibenkan, ilmend n twila.

Talɣut yettwaseknen

Sin n yibenkan

Tafelwit i sin n yibenkan deg tagrumma Rondo.

Timeẓliyin yezdin tagrumma Panvisus Rondo. Tamsisɣelt n tfelwit, asenqed n teɣzi d yimerna ilaq ad ttwasentmen.

www.fiso.info · 2026-10-18

Tafelwit ilmend n ufran

Tafelwit s talɣa n V neɣ s uselled, ilmend n temsisɣelt.

www.fiso.info · 2026-10-18

Teɣzi ilmend n temsisɣelt

Tafelwit s teɣzi tusbitt neɣ s useggem s trisiti, ilmend n lqem.

www.fiso.info · 2026-10-18

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

Huvitz | ASEM

Réfracteur HDR-7100P



Arifraktur i ukayad asaguran seg lbeɛd d lqerb.

Talɣut yettwaseknen

Rifraksyun tasagurant

SPH seg −29 ar +26,75 D i ukayad n tmugli.

huvitz.com · 2026-10-05

Tamugli n lqerb

Lbeɛd n umahil yettuseggam seg 35 ar 70 cm.

huvitz.com · 2026-10-05

Aseɣti aprizmatik

Tazmert taprizmatikt alamma 20 Δ.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Haag-Streit | ASEM

Lampe à fente BM-900



Taftilt s tecraḍ, amikruskuḍ Kepler d yisemeɣriyen 10×/16×.

Talɣut yettwaseknen

Taḍwant n umikruskuḍ Kepler

Sin n yisemeɣriyen 10× d 16× i tmugli s tecraḍ.

haag-streit.com · 2026-10-17

Tafat LED n tecraḍ

Tafat LED i tmugli s teftilt s tecraḍ.

Lqem LED

jp.haag-streit.com · 2026-10-17

Asers s yiwen n ufus

Asers n 3 n yigellusen n yibenk s uɣil n yiwen n ufus.

jp.haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Huvitz | ASEM

Frontofocometre Digital



Akat n lentiyat n tira d unermis s tetiknuluji Hartmann.

Talɣut yettwaseknen

Tasleḍt Hartmann

Akat n SPH, CYL d tzemmart taprizmatikt.

huvitz.com · 2026-10-05

Lentiyat n unermis

Akat n lentiyat n unermis iquranen d yilwiyen.

huvitz.com · 2026-10-05

Azilal aḍwan wessiɛen

SPH seg −25 ar +25 D; CYL alamma 10 D.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

HIIS-1 (HIS-1)



Aseɣẓan i twali n yikayaden i d-yekkan seg yibenkan imṣadan.

Talɣut yettwaseknen

Tamuɣli n yiwen ukayad

Askan n yiwen ukayad i twali n yisefka-ines.

Tawriqt tessemras ismawen HIS-1 d HIIS-1. Lqem, turagin, izegrar d yibenk nniḍen yeqqnen ilaq ad ttwasentmen.

optomedical.eu · 2026-10-18

Tamuɣli n sdis n yikayaden

Tamuɣli n sdis n yikayaden i userwes n yigmaḍ.

optomedical.eu · 2026-10-18

Ambaddal s DICOM

Aneggaf DICOM i umbaddal akked yinagrawen imṣadan.

optomedical.eu · 2026-10-18

Wali talqayt

Wali taɣbalut n umfaras · optomedical.eu

Haag-Streit | ASEM

AT 900



Amektil n uẓeẓẓu n tiṭ s umendad n Goldmann n usbeṭṭeḥ, seg tagrumma AT 900.

Talɣut yettwaseknen

Amendad n Goldmann

Aktil n uẓeẓẓu daxel n tiṭ s usbeṭṭeḥ.

Timeẓliyin yezdin tagrumma AT 900. Aseɣwen R, T neɣ BQ d umṣada akked tafat n tfenda ilaq ad ttwasentmen.

uk.haag-streit.com · 2026-10-18

Azilal n uktil

Azilal n uẓeẓẓu iɣef yella yisemli: seg 0 ar 80 mmHg.

uk.haag-streit.com · 2026-10-18

Aseɣwen i ilaq ad yettwasentem

Aseɣwen R, T neɣ BQ d umṣada ilaq ad ttwasentmen.

uk.haag-streit.com · 2026-10-18

Wali talqayt

Wali taɣbalut n umfaras · uk.haag-streit.com

Huvitz | ASEM

Auto-réfkeratomètre HRK-1



Rifraksyun, keratumetri d tallalt i usmṣada n lentiyat n unermis.

Talɣut yettwaseknen

Rifraksyun d lqerniya

Rifraktumetri d keratumetri deg yiwen n yibenk.

huvitz.com · 2026-10-05

Asmṣada n lentiyat

Tallalt i usmṣada n lentiyat n unermis.

huvitz.com · 2026-10-05

Aḍfar ɣef teɣzi

Aḍfar n umḍiq ɣef teɣzi ɣef ±15 mm.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · www.huvitz.com

Huvitz | ASEM

Frontofocometre HLM-7000P



Akat n lentiyat n tira d lentiyat n unermis.

Talɣut yettwaseknen

Akat n lentiyat n tira

Akat n SPH, CYL d tzemmart taprizmatikt.

huvitz.com · 2026-10-05

Lentiyat n unermis

Akat n lentiyat n unermis iquranen d yilwiyen.

huvitz.com · 2026-10-05

Azilal aḍwan wessiɛen

SPH seg −25 ar +25 D; CYL alamma 10 D.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

Projecteur de test HCP-7000



Asnay LED i yikayaden n tmugli deg uswir.

Talɣut yettwaseknen

Ikayaden n tmugli yemgaraden

41 n yizmal n tmugli d 34 n tgelmamin n uskan.

huvitz.com · 2026-10-05

Asnay LED

Aɣbalu n tafat LED i usnay n yizmal n tmugli.

huvitz.com · 2026-10-05

Lbeɛd n ukayad

Asnay i yilbɛad n ukayad seg 2,5 ar 8 m.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

Ecran Numérique HDC-9100N/PF



Askan umḍin n yizmal n tmugli ɣef ugdil LCD Full HD.

Talɣut yettwaseknen

Izmal n tmugli Full HD

Agdil LCD n 24″, tafrayt 1920 × 1080.

huvitz.com · 2026-10-05

Lbeɛd n ukayad

Askan n yikayaden n tmugli ɣer lbeɛd n 6 m.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Champ visuel HVF-100



Akat n wurti n tmugli s perimetri ur nettbeddil ara, amellal ɣef umellal.

Talɣut yettwaseknen

Perimetri ur nettbeddil ara

Akayad n wurti n tmugli, amellal ɣef umellal.

huvitz.com · 2026-10-05

Urti alamma 90°

Akayad aperimetrik n wurti alamma 90°.

huvitz.com · 2026-10-05

Amattay Goldmann III

Tiddi n umattay Goldmann III i yikayaden n wurti n tmugli.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

Tonomètre à air HNT-1/1P



Tonumetri war anermis i ukat n tḍeɣṭ deg tiṭ.

Talɣut yettwaseknen

Tonumetri war anermis

Akat n tḍeɣṭ deg tiṭ war anermis.

huvitz.com · 2026-10-05

Azilal n tḍeɣṭ

Akat n tḍeɣṭ deg uzilal n 0–60 mmHg.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

Huvitz | ASEM

Topographe Optique Cornéen Complet



Topografi n lqerniya s keratumetri d tesleḍt Zernike.

Talɣut yettwaseknen

Takarḍa n lqerniya

Takarḍa tatopugrafit n tmaẓirt n lqerniya.

huvitz.com · 2026-10-05

Tasleḍt Zernike

Tasleḍt n isembeddilen n lqerniya s wawalen Zernike.

huvitz.com · 2026-10-05

Keratumetri deg yiwen

Akat n ukrif n lqerniya deg yiwen n ukayad.

huvitz.com · 2026-10-05

Wali talqayt

Wali taɣbalut n umfaras · huvitz.com

F.I.S.O. | ASEM

Unité 2 postes



Nermes ASEM i yisallen yettwasenqden ɣef temselɣut-a.

Tulmisin titiknikanin ilaq ad ttwasentmen.

Wali talqayt

Wali taɣbalut n umfaras · www.fiso.info

Haag-Streit | ASEM

EIBOS 2



Aferdis n tmugli wessiɛen i umikruskuḍ n ugzay yemṣadan.

Talɣut yettwaseknen

Tamugli wessiɛen war anermis

Tamugli wessiɛen n retina s umikruskuḍ n ugzay.

Amikruskuḍ yemṣadan ilaq

haag-streit.com · 2026-10-17

Amuttay n tugna deg yiwen

Alemmes d umuttay n tugna deg yiwen n unagraw.

haag-streit.com · 2026-10-17

Alemmes s ufus

Iɣallen n tama i useggem n ulemmes s ufus.

haag-streit.com · 2026-10-17

Wali talqayt

Wali taɣbalut n umfaras · haag-streit.com

Huvitz HS-7000 — taftilt s tecraḍ



Taftilt s tecraḍ tagalilit d semmus n yisemeɣriyen.

Talɣut yettwaseknen

Semmus n yisemeɣriyen

6×, 10×, 16×, 25×, 40× i usmṣada n tmugli.

huvitz.com · 2026-10-06

Tacraḍ n ukayad

Aldi n tecraḍ alamma 14 × 14 mm.

huvitz.com · 2026-10-06

Tistayin n tmugli

Tastayt tazegzawt cobalt d tazegzawt i ukayad.

huvitz.com · 2026-10-06

Wali talqayt

Wali taɣbalut n umfaras · www.huvitz.com

Huvitz | ASEM

Huvitz HS-5000 — taftilt s tecraḍ



Taftilt s tecraḍ d tafat n tburt, s tḍwant tagalilit.

Talɣut yettwaseknen

Semmus n yisemeɣriyen

6,5×, 11×, 16×, 23,2×, 39,3× i usmṣada n tmugli.

huvitz.com · 2026-10-06

Tacraḍ n ukayad

Aldi n tecraḍ alamma 14 × 14 mm.

huvitz.com · 2026-10-06

Tafat n tburt

Aḍran n tafat s yiswiren seg 0° ar 20°.

huvitz.com · 2026-10-06

Wali talqayt

Wali taɣbalut n umfaras · www.huvitz.com

Huvitz HS-5500 — taftilt s tecraḍ



Tafat deg yiwen d tḍwant tagalilit i teftilt s tecraḍ.

Talɣut yettwaseknen

Semmus n yisemeɣriyen

6,5×, 11×, 16×, 23,2×, 39,3× i usmṣada n tmugli.

huvitz.com · 2026-10-06

Tacraḍ n ukayad

Aldi n tecraḍ alamma 14 × 14 mm.

huvitz.com · 2026-10-06

Tafat deg yiwen

Aḍran n tafat ur nettbeddil ara s yiswiren seg 0° ar 20°.

huvitz.com · 2026-10-06

Wali talqayt

Wali taɣbalut n umfaras · www.huvitz.com

Huvitz HDC-7000 — anagraw n yikayaden n tmugli



Ikayaden umḍinen i rifraksyun d tmugli n snat n wallen.

Talɣut yettwaseknen

Tamugli n snat n wallen

Ikayaden n stereopsis, asdukel d umgired n usbadu.

www.huvitz.com · 2026-10-06

Ikayaden yettusmṣadan

Amgired, tigelmamin, azeggaɣ/azegzaw d uskan s lemri.

www.huvitz.com · 2026-10-06

Asbeddi ɣef PC

Ikayaden ɣef PC d ugdil yemmṣadan i ilaq ad ttwasbeggnen.

PC d ugdil berra

www.huvitz.com · 2026-10-06

Wali talqayt

Wali taɣbalut n umfaras · www.huvitz.com

Huvitz CIT-4000 — tafelwit s umutur



Tafelwit s umutur i usebded n yibenk n ukayad.

Talɣut yettwaseknen

Aseggem s umutur

Aseggem elektrik n teɣzi n yibenk.

www.huvitz.com · 2026-10-06

Amussu deg uswir

Tirwiḍin i usenfel n tfelwit gar yimukan n ukayad.

www.huvitz.com · 2026-10-06

Asenqed yettusmṣadan

Asenqed s ufus neɣ s uḍar ilmend n lqem.

Ilmend n twila

www.huvitz.com · 2026-10-06

Wali talqayt

Wali taɣbalut n umfaras · www.huvitz.com