Tayunin n usekyed n tɣuri

Easyvisus RM64 — 2 n yibenk

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Easyvisus RM64 i sin n yibenkan yesdukkel ṭṭabla s teɣzi s trisiti d lbiru s tjabbuḍt i tfelwit n 236 n tlentiyin n ukayad. F.I.S.O. yeskan leqwamat yettuzzuyun, s neɣ war tajabbuḍt tamernant, d ṭṭwabel yettselliden s ufus neɣ s trisiti. Agdil LCD d ujuyestik ḥekkmen deg users n ṭṭabla d ukersi PRIMA. Tamsisɣelt s timmad-is, ibenkan d yimallalen ilaq ad ttwasegzin deg tesɣunt ASEM.



Easyvisus RM64 · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-09

> Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Asebded elektrik

Aseggem elektrik n tfelwit d ukersi.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Asewḥel elektromagnitik

Asewḥel n tfelwit d uɣil n urifraktur.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Suddes amkan Easyvisus RM64 — 2 n yibenkan

Tugna n umeskar teskan tawila s timmad-is d unekcum ɣer yibenkan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30


Serres tafelwit n umahil

Amussu n tfelwit iɣef yella yisemli yeshel abeddel seg yibenk ɣer wayeḍ deg Easyvisus RM64 — 2 n yibenkan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Heggi asbeddi

Azenziɣ unsib yeskan aḥric ilaqen deg wadeg, yettɛawan aheggi n usbeddi n Easyvisus RM64 — 2 n yibenkan.



Easyvisus RM64 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Ergonomi

Asbeddi yettwaxedmen ilmend n umkan n ukayad.

Ibenkan yettwasuddesen

Ibenkan ttaren amkan-nsen deg yiwet n tayunt tusligt.

Iseɣtayen n umkan

Imussuten yellan seshlen asgunfu n twila.

Tiwuriwin deg usnas

Tayunt n rifraksyun

Iɣbula d yisemliyen

I wakken ad yishil ubrid n umahil

Ṭṭabla iwulmen i yibenkan yettwafrenen

400600 tesɛa ṭṭabla yettuzzuyun war tajabbuḍt n ṭṭabla; 400630 ternu tajabbuḍt s tfelwit n usɣar i 148 n tlentiyin n ukayad. 400650 tesɛa aselled s ufus, 400655 aselled s trisiti. Akk sɛant aseggem n teɣzi s trisiti d lbiru s tjabbuḍt i tfelwit n 236 n tlentiyin. 400630 tesɛa ihi snat n tjabbuḍin yemgaraden. Tilentiyin d yibenkan ur ttwaḥsaben ara dduklen.

Iɣbula d yisemliyen

Inadiyen d users s yisemliyen

Agdil LCD yeğğa aseggem n lugu. Ajuyestik yettḥakkam deg usali d usrus n ukersi d ṭṭabla d usekkeṛ. Ameskar yeskan asekkeṛ iliktrumagnitik n ṭṭabla yettuzzuyun d ufus n ṭṭabla. Amutur n teɣzi yemgarad ɣef tuzzya neɣ aselled yettwafren i temsisɣelt.

Iɣbula d yisemliyen

I useqdec n yal ass

Akersi amagnu d yimallalen n ufran

PRIMA amagnu yesɛa aseggem n teɣzi s trisiti d taduli skay taberkant, war asersi n yiḍarren. Asersi n yiḍarren i yettalayen, tidulin-nniḍen d ubeddel s wanaw-nniḍen n ukersi d ifranen iɣef llan yisemliyen. Leqwamat n ukersi ilaq ad ttwaseɣten deg tesɣunt ASEM. Tawlaft n tazwara tezmer ad tesken asersi n yiḍarren, afus n urifraktyur neɣ imallalen-nniḍen n ufran; askan-nsen ur iḍemmen ara dduklen. Ibenkan d yimallalen ilaq ad ttwasegzin s yinagrawen yemgaraden.

Iɣbula d yisemliyen

Tafawt d trisiti i yibenkan

Tagejdit n laluminyum tesɛa tafawt LED n twennaḍt i yettuseggamen, taftilt LED n tɣuri d tfelwit n users i usenfal n yikayaden. Tayunt tagejdant tesɛa trisiti n utransfurmatur n taftilt n tfexsiḍt. Llsas n wuzzal yetteḍlan yesɛa iḍarren i yettuseggamen. Iɣbula n trisiti imernan i yibenkan llan deg yimallalen: ur ttwaḥsaben ara d ayen yedduklen s wudem awurman. Ulac tazwert n trisiti neɣ asnagar n tayunt meṛṛa i d-yettwakkasen seg leḥsab.

Iɣbula d yisemliyen



Easyvisus RM64 · Tugna

Iɣbula d yisemliyen



Easyvisus RM64 · Tugna

Iɣbula d yisemliyen

Asbeddi d ufran n twila

Lqem s uselled s trisiti 400655 yella i yefus kan. Twiliwin-nniḍen zemrent ad ttwafranent i yizelmaḍen: 400600S, 400630S d 400650S. Izenziɣen semgaraden ṭṭwabel yettuzzuyun d tid yettselliden; ikaten-nsent ur skanen ara yiwet n tsekta tamatut yemṣadan. Ulac taẓayt, asnagar n trisiti neɣ tilisa n teɣzi i d-yettwakkasen seg yizenziɣen-a. Asbeddi d temsisɣelt ilaq ad ttwaseɣten akked ASEM.

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
F.I.S.O.

Anaw
Easyvisus RM64 — 2 n yibenk

Tatiknuluji
Tayunt n rifraksyun

Aseqdec akliniki
Itiknulujiwin tusligin

Amḍan n yibenkan
2

Iɣbula d yisemliyen

Teɣzi n ukersi
Aseggem s trisiti; akersi Prima amagnu

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Teɣzi n ṭṭabla
Aseggem s trisiti

Iɣbula d yisemliyen

Asekkeṛ n ṭṭabla
Ṭṭabla d ufus: asekkeṛ aliktrumagnitik

Iɣbula d yisemliyen

Amussu n ṭṭabla
Tuzzya neɣ amussu usrid ilmend n temsisɣelt

Afran n twila n umeskar; ulac afran awurman i yettwaḥsaben yettwasebded.

Iɣbula d yisemliyen

Tafawt n umkan n umahil

LED n twennaḍt · LED n tɣuri

Iɣbula d yisemliyen

Inadiyen

Agdil LCD; ajuyestik i useggem d usekkeṛ

Iɣbula d yisemliyen


Llsas

Uzzal yetteḍlan; iḍarren yettuseggamen

Iɣbula d yisemliyen

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Afran n twila n umeskar; ulac afran awurman i yettwaḥsaben yettwasebded.



Easyvisus RM64 · Tugna

Iɣbula d yisemliyen

Asider d yiɣbula

Isemliyen n umeskar

PDF · 15,672 kB

Iɣbula d yisemliyen

- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-unit-2-instruments_451/
  Aɣbalu yettwasentem ass n 2026-10-05
- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-unit-2-instruments_451/
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-unit-2-instruments_451/
  Aɣbalu yettwasentem ass n 2026-10-17

- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-unit-2-instruments_451/
  Aɣbalu yettwasentem ass n 2026-08-30
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef Easyvisus RM64 — 2 n yibenk

D acu Amḍan n yibenkan iɣef llan yisemliyen n Easyvisus RM64 — 2 n yibenk?

2

Iɣbula d yisemliyen

D acu Teɣzi n ukersi iɣef llan yisemliyen n Easyvisus RM64 — 2 n yibenk?

Aseggem s trisiti; akersi Prima amagnu — Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anta tawila i ilaq ad sentmeɣ i Easyvisus RM64 — 2 n yibenk?

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anda ara afeɣ isemliyen n Easyvisus RM64 — 2 n yibenk?

Isemliyen n umeskar (PDF) — https://asem-medical.com/storage/catalogue-enrichment/fiso-20261006/fiso-original-2025.pdf

Iɣbula d yisemliyen

Anwa imussuten n tfelwit i yettwasegzan?

RM64 n sin n yibenkan yella s tfelwit yettuzuzuren 400600 neɣ 400630, azleg s ufus 400650, neɣ azleg atrisiti 400655. Akk ttwasfuklen s useggem atrisiti n teɣzi.

Iɣbula d yisemliyen

Azleg atrisiti yella i yizelmaḍen?

Tawriqt teḥbes 400655 atrisiti deg lqem ayeffus kan. Tiwila 400600, 400630 d 400650 zemrent ad sɛunt lqem azelmaḍ d afran.

Iɣbula d yisemliyen

Aḥraz kifkif deg yal lqem?

Tansa n umahil tettwahegga i tfelwit n 236 n lentes n usekyed. 400630 yernu tajabbuḍt ɣer tfelwit n tayunt i 148 n lentes; yesɛa ihi snat n tejabbuḍin. Ikitaten n lentes ilaq ad ttusegzin iman-nsen.

Iɣbula d yisemliyen

Ifarisen icudden



Panvisus Rondo

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt Panvisus Rondo i sin n yibenkan, ilmend n twila.

Wali afaris



Easyvisus Curve

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 2 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 3 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 3 n yibenk n ukayad.

Wali afaris

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Tetiknulujiyin tusliḍin

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: d87fccce45c8