Tayunin n usekyed n tɣuri

Easyvisus RM64 Compact — 2 n yibenk

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Easyvisus RM64 Compact i sin n yibenkan d tayunt s ṭṭabla n teɣzi i yettuseggamen s trisiti. Tawsit tesɛa tuzzya, aselled s ufus neɣ s trisiti, d tjabbuḍt n ṭṭabla deg yiwet n twila tusligt. Tesɛa agdil LCD, ajuyestik d ukersi PRIMA s trisiti, war lbiru yemgaraden i tfelwit n 236 n tlentiyin n RM64 ur nelli d Compact. Tamsisɣelt iwulmen i tklinik ilaq ad tettwaseɣti akked ASEM.



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-09

Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Anermis s ukersi s tirwiḍin

Tafelwit tettildi i ukayad n umuḍin ɣef ukersi s tirwiḍin.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Suddes amkan Easyvisus RM64 Compact — 2 n yibenkan

Tugna n umeskar teskan tawila s timmad-is d unekcum ɣer yibenkan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30


Serres tafelwit n umahil

Amussu n tfelwit iɣef yella yisemli yeshel abeddel seg yibenk ɣer wayeḍ deg Easyvisus RM64 Compact — 2 n yibenkan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Heggi asbeddi

Azenziɣ unsib yeskan aḥric ilaqen deg wadeg, yettɛawan aheggi n usbeddi n Easyvisus RM64 Compact — 2 n yibenkan.



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Ergonomi

Asbeddi yettwaxedmen ilmend n umkan n ukayad.

Ibenkan yettwasuddesen

Ibenkan ttaren amkan-nsen deg yiwet n tayunt tusligt.

Iseɣtayen n umkan

Imussuten yellan seshlen asgunfu n twila.

Tiwuriwin deg usnas

Tayunt n rifraksyun

Iɣbula d yisemliyen

I wakken ad yishil ubrid n umahil

Ṭṭabla iwulmen i yibenkan yettwafrenen

400900 d 400930 ttzuzzuyent; 400930 kan i yesɛan tajabbuḍt s tfelwit n usɣar i 148 n tlentiyin n ukayad. 400950 tesɛa aselled s ufus, 400955 s trisiti. Akk sɛant teɣzi n ṭṭabla s trisiti. Ameskar yettak adday n ṭṭabla n 76 cm d uldi i ukayad n umdan s ukersi s rrwaḍi. Tilentiyin ilaq ad ttwasegzin, adday-a ur yeskan ara teɣzi tamatut n tayunt.

Iɣbula d yisemliyen

Inadiyen d users s yisemliyen

Agdil LCD yeğğa aseggem n lugu. Ajuyestik yettḥakkam deg usali d usrus n ukersi d ṭṭabla d usekkeṛ. Ameskar yeskan asekkeṛ iliktrumagnitik n ṭṭabla yettuzzuyun d ufus n ṭṭabla. Amutur n teɣzi yemgarad ɣef tuzzya neɣ aselled yettwafren i temsisɣelt.

Iɣbula d yisemliyen

I useqdec n yal ass

Akersi amagnu d yimallalen n ufran

PRIMA amagnu yesɛa aseggem n teɣzi s trisiti d taduli skay taberkant, war asersi n yiḍarren. Asersi n yiḍarren i yettalayen, tidulin-nniḍen d ubeddel s wanaw-nniḍen n ukersi d ifranen iɣef llan yisemliyen. Leqwamat n ukersi ilaq ad ttwaseɣten deg tesɣunt ASEM. Tawlaft n tazwara tezmer ad tesken asersi n yiḍarren, afus n urifraktyur neɣ imallalen-nniḍen n ufran; askan-nsen ur iḍemmen ara dduklen. Ibenkan d yimallalen ilaq ad ttwasegzin s yinagrawen yemgaraden.

Iɣbula d yisemliyen

Tafawt d trisiti i yibenkan

Tagejdit n laluminyum tesɛa tafawt LED n twennaḍt i yettuseggamen, taftilt LED n tɣuri d tfelwit n users i usenfal n yikayaden. Tayunt tagejdant tesɛa trisiti n utransfurmatur n taftilt n tfexsiḍt. Llsas n wuzzal yetteḍlan yesɛa iḍarren i yettuseggamen. Iɣbula n trisiti imernan i yibenkan llan deg yimallalen: ur ttwaḥsaben ara d ayen yedduklen s wudem awurman. Ulac tazwert n trisiti neɣ asnagar n tayunt meṛṛa i d-yettwakkasen seg leḥsab.

Iɣbula d yisemliyen



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Asbeddi d ufran n twila

Aselled s trisiti 400955 yettunefk i yefus kan. Twiliwin-nniḍen sɛant leqwamat n ufran i yizelmaḍen 400900S, 400930S d 400950S. Lbiru yemgaraden s tjabbuḍt i 236 n tlentiyin n RM64 ur nelli d Compact maččі d ibenk amagnu n lqem-a. Talɣiwin yettuzzuyun d tid yettselliden yemgaradent. Ulac taẓayt, asnagar neɣ teɣzi tafellayt i yettwasnulfan; adeg n umussu ilaq ad yettwaseɣti akked ASEM.

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
F.I.S.O.

Anaw
Easyvisus RM64 Compact — 2 n yibenk

Tatiknuluji
Tayunt n rifraksyun

Aseqdec akliniki
Itiknulujiwin tusligin

Amḍan n yibenkan
2

Iɣbula d yisemliyen

Teɣzi n ukersi
Aseggem s trisiti; akersi Prima amagnu

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Teɣzi n ṭṭabla
Aseggem s trisiti

Iɣbula d yisemliyen

Asekkeṛ n ṭṭabla
Ṭṭabla d ufus: asekkeṛ aliktrumagnitik

Iɣbula d yisemliyen

Amussu n ṭṭabla
Tuzzya neɣ amussu usrid ilmend n temsisɣelt

Afran n twila n umeskar; ulac afran awurman i yettwaḥsaben yettwasebded.

Iɣbula d yisemliyen

Tafawt n umkan n umahil

LED n twennaḍt · LED n tɣuri

Iɣbula d yisemliyen

Inadiyen

Agdil LCD; ajuyestik i useggem d usekkeṛ

Iɣbula d yisemliyen


Llsas

Uzzal yetteḍlan; iḍarren yettuseggamen

Iɣbula d yisemliyen

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Afran n twila n umeskar; ulac afran awurman i yettwaḥsaben yettwasebded.



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Asider d yiɣbula

Isemliyen n umeskar

PDF · 15,672 kB

Iɣbula d yisemliyen

- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-2-instruments_496/
  Aɣbalu yettwasentem ass n 2026-10-05
- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-2-instruments_496/
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-2-instruments_496/
  Aɣbalu yettwasentem ass n 2026-10-17

- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-2-instruments_496/
  Aɣbalu yettwasentem ass n 2026-08-30
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef Easyvisus RM64 Compact — 2 n yibenk

D acu Amḍan n yibenkan iɣef llan yisemliyen n Easyvisus RM64 Compact — 2 n yibenk?

2

Iɣbula d yisemliyen

D acu Teɣzi n ukersi iɣef llan yisemliyen n Easyvisus RM64 Compact — 2 n yibenk?

Aseggem s trisiti; akersi Prima amagnu — Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anta tawila i ilaq ad sentmeɣ i Easyvisus RM64 Compact — 2 n yibenk?

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anda ara afeɣ isemliyen n Easyvisus RM64 Compact — 2 n yibenk?

Isemliyen n umeskar (PDF) — https://asem-medical.com/storage/catalogue-enrichment/fiso-20261006/fiso-original-2025.pdf

Iɣbula d yisemliyen

Anwa afran n tfelwit i d-yettak Compact n sin n yibenkan?

Tawriqt tefrez tifelwin yettuzuzuren 400900 d 400930, azleg s ufus 400950 d uzleg atrisiti 400955. 400930 yesɛa tajabbuḍt i tfelwit n 148 n lentes n usekyed.

Iɣbula d yisemliyen

Anwa tilisa i yerzan lqem azelmaḍ?

400955 s uzleg atrisiti yettwasegza s talɣa tayeffust kan. 400900, 400930 d 400950 zemren ad d-ttunefken s lqem azelmaḍ d afran.

Iɣbula d yisemliyen

Tawriqt teskan anekcum s ukersi s ṛṛwāḍ?

Ih, teskan alday n tfelwit, s teɣzi taddayt n 76 cm, i umuḍin i yeqqimen deg ukersi-ines s ṛṛwāḍ. Tazmert-a ilaq ad tettwasenqed deg usbeddi yettwafernen.

Iɣbula d yisemliyen

Ifarisen icudden



Panvisus Rondo

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt Panvisus Rondo i sin n yibenkan, ilmend n twila.

Wali afaris



Easyvisus Curve

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 2 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 3 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 3 n yibenk n ukayad.

Wali afaris

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Tetiknulujiyin tusliḍin

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: c64d631e6a4f