Tayunin n usekyed n tɣuri

Easyvisus RM64 Compact — 3 n yibenk

Tayunt n rifraksyun i 3 n yibenk n ukayad.

Easyvisus RM64 Compact i kraḍ n yibenkan yesɛa ṭṭabla s teɣzi s trisiti, tuzzya s ufus neɣ asers awurman ilmend n temsisɣelt. S umgired ɣef RM64 s lbiru, Compact-a ur yeskan ara lbiru yemgaraden s tjabbuḍt i tfelwit n 236 n tlentiyin. Agdil LCD, ajuyestik d ukersi PRIMA s trisiti ssudusen akayad. Ifranen d yibenkan ilaq ad ttwasegzin akked ASEM.



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-09

Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Amkan i 3 n yibenk

Tafelwit tettusebda i 3 n yibenk n ukayad.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Anermis s ukersi s tirwiḍin

Tafelwit tettildi i ukayad n umuḍin ɣef ukersi s tirwiḍin.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Suddes amkan Easyvisus RM64 Compact — 3 n yibenkan

Tugna n umeskar teskan tawila s timmad-is d unekcum ɣer yibenkan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30


Serres tafelwit n umahil

Amussu n tfelwit iɣef yella yisemli yeshel abeddel seg yibenk ɣer wayeḍ deg Easyvisus RM64 Compact — 3 n yibenkan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Heggi asbeddi

Azenziɣ unsib yeskan aḥric ilaqen deg wadeg, yettɛawan aheggi n usbeddi n Easyvisus RM64 Compact — 3 n yibenkan.



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Ergonomi

Asbeddi yettwaxedmen ilmend n umkan n ukayad.

Ibenkan yettwasuddesen

Ibenkan ttaren amkan-nsen deg yiwet n tayunt tusligt.

Iseɣtayen n umkan

Imussuten yellan seshlen asgunfu n twila.

Tiwuriwin deg usnas

Tayunt n rifraksyun

Iɣbula d yisemliyen

I wakken ad yishil ubrid n umahil

Ṭṭabla iwulmen i yibenkan yettwafrenen

400980 tesɛa ṭṭabla s tuzzya s ufus i kraḍ n yibenkan, 400985 s users awurman. Sin n leqwamat sɛan aseggem n teɣzi s trisiti. Ṭṭabla tezmer ad tettwaldi s lekmal i ukayad n umdan yellan deg ukersi s rrwaḍi; teɣzi taddayt iɣef yella yisemli d 76 cm. Adday-a mačči d teɣzi tamatut n tayunt neɣ d tilisa meṛṛa n useggem.

Iɣbula d yisemliyen

Inadiyen d users s yisemliyen

Agdil LCD yeğğa aseggem n lugu. Ajuyestik yettḥakkam deg usali d usrus n ukersi d ṭṭabla d usekkeṛ. Ameskar yeskan asekkeṛ iliktrumagnitik n ṭṭabla yettuzzuyun d ufus n ṭṭabla. Amutur n teɣzi yemgarad ɣef tuzzya neɣ aselled yettwafren i temsisɣelt.

Iɣbula d yisemliyen

I useqdec n yal ass

Akersi amagnu d yimallalen n ufran

PRIMA amagnu yesɛa aseggem n teɣzi s trisiti d taduli skay taberkant, war asersi n yiḍarren. Asersi n yiḍarren i yettalayen, tidulin-nniḍen d ubeddel s wanaw-nniḍen n ukersi d ifranen iɣef llan yisemliyen. Leqwamat n ukersi ilaq ad ttwaseɣten deg tesɣunt ASEM. Tawlaft n tazwara tezmer ad tesken asersi n yiḍarren, afus n urifraktyur neɣ imallalen-nniḍen n ufran; askan-nsen ur iḍemmen ara dduklen. Ibenkan d yimallalen ilaq ad ttwasegzin s yinagrawen yemgaraden.

Iɣbula d yisemliyen

Tafawt d trisiti i yibenkan

Tagejdit n laluminyum tesɛa tafawt LED n twennaḍt i yettuseggamen, taftilt LED n tɣuri d tfelwit n users i usenfal n yikayaden. Tayunt tagejdant tesɛa trisiti n utransfurmatur n taftilt n tfexsiḍt. Llsas n wuzzal yetteḍlan yesɛa iḍarren i yettuseggamen. Iɣbula n trisiti imernan i yibenkan llan deg yimallalen: ur ttwaḥsaben ara d ayen yedduklen s wudem awurman. Ulac tazwert n trisiti neɣ asnagar n tayunt meṛṛa i d-yettwakkasen seg leḥsab.

Iɣbula d yisemliyen



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Asbeddi d ufran n twila

Lqem n ufus yezmer ad yettwafren i yizelmaḍen s temsisɣelt 400980S. 400985 s users awurman yella i yefus kan. Izenziɣen qqnen ɣer twiliwin-a yemgaraden; ulac yiwet n tsekta tamatut i d-yettwakkasen. Lbiru s tjabbuḍt i 236 n tlentiyin n RM64 ur nelli d Compact ur yeddukel ara deg yibenkan imagnuten i d-yettwaskanen dagi. Taẓayt, asnagar d teɣzi tafellayt qqimen war isemliyen deg yisebtar-a.

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
F.I.S.O.

Anaw
Easyvisus RM64 Compact — 3 n yibenk

Tatiknuluji
Tayunt n rifraksyun

Aseqdec akliniki
Itiknulujiwin tusligin

Amḍan n yibenkan
3

Iɣbula d yisemliyen

Teɣzi n ukersi
Aseggem s trisiti; akersi Prima amagnu

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Teɣzi n ṭṭabla
Aseggem s trisiti

Iɣbula d yisemliyen

Asekkeṛ n ṭṭabla
Ṭṭabla d ufus: asekkeṛ aliktrumagnitik

Iɣbula d yisemliyen

Amussu n ṭṭabla
Tuzzya s ufus; asers awurman ilmend n temsisɣelt

Afran n twila n umeskar; ulac afran awurman i yettwaḥsaben yettwasebded.

Iɣbula d yisemliyen

Tafawt n umkan n umahil

LED n twennaḍt · LED n tɣuri

Iɣbula d yisemliyen

Inadiyen

Agdil LCD; ajuyestik i useggem d usekkeṛ

Iɣbula d yisemliyen


Llsas

Uzzal yetteḍlan; iḍarren yettuseggamen

Iɣbula d yisemliyen

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Afran n twila n umeskar; ulac afran awurman i yettwaḥsaben yettwasebded.



Easyvisus RM64 Compact · Tugna

Iɣbula d yisemliyen

Asider d yiɣbula

Isemliyen n umeskar

PDF · 15,672 kB

Iɣbula d yisemliyen

- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-3-instruments_499/
  Aɣbalu yettwasentem ass n 2026-10-05
- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-3-instruments_499/
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-3-instruments_499/
  Aɣbalu yettwasentem ass n 2026-10-17

- https://www.fiso.info/en/products/refraction_units/easyvisus-rm64-compact-3-instruments_499/
  Aɣbalu yettwasentem ass n 2026-08-30
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef Easyvisus RM64 Compact — 3 n yibenk

D acu Amḍan n yibenkan iɣef llan yisemliyen n Easyvisus RM64 Compact — 3 n yibenk?

3

Iɣbula d yisemliyen

D acu Teɣzi n ukersi iɣef llan yisemliyen n Easyvisus RM64 Compact — 3 n yibenk?

Aseggem s trisiti; akersi Prima amagnu — Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anta tawila i ilaq ad sentmeɣ i Easyvisus RM64 Compact — 3 n yibenk?

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anda ara afeɣ isemliyen n Easyvisus RM64 Compact — 3 n yibenk?

Isemliyen n umeskar (PDF) — https://asem-medical.com/storage/catalogue-enrichment/fiso-20261006/fiso-original-2025.pdf

Iɣbula d yisemliyen

Anwa anekcum i yettwasegzan i umuḍin s ukersi s ṛṛwāḍ?

Tawriqt teskan tafelwit i yeldin s teɣzi taddayt n 76 cm i ukayad n umuḍin i yeqqimen deg ukersi-ines s ṛṛwāḍ. Asbeddi d twila yettwafernen ilaq ad ttwasenqden akked ASEM.

Iɣbula d yisemliyen

Anwa tuzzyiwin i d-yettunefken?

400980 yessemras tuzzya s ufus; 400985 yettak aseggem awurman. Sin yid-sen rzan tafelwit n kraḍ n yibenkan s teɣzi atrisitant.

Iɣbula d yisemliyen

Lqem azelmaḍ yella deg sin n yinawen?

Ala. Lqem azelmaḍ 400980S yettwasegza i tuzzya s ufus. Aseggem awurman yettwasegza deg lqem ayeffus kan.

Iɣbula d yisemliyen

Ifarisen icudden



Panvisus Rondo

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt Panvisus Rondo i sin n yibenkan, ilmend n twila.

Wali afaris



Easyvisus Curve

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 2 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 3 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 3 n yibenk n ukayad.

Wali afaris

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Tetiknulujiyin tusliḍin

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: 85f19a3022ea