Tayunin n usekyed n tɣuri

Easyvisus Two

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Easyvisus Two d tayunt n rifraksyun s tfelwit V n teɣzi ur nettbeddil ara i sin n yibenkan. Amekanizm uslig yeslemmas ibenkan, afus n ṭṭabla yettusekkeṛ s umekanizm. Akersi PRIMA s teɣzi s trisiti, taftilt LED n tɣuri d tjabbuḍt i yettwakkasen, yettḥadar ɣef uɣebbar, d tulmisin timagnutin s yisemliyen. Tafelwit i wer nettbeddil ara ur ilaq ara ad tettwammel s umutur deg ibeddi.



Easyvisus Two · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-09

> Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Amkan i 2 n yibenk

Tafelwit tettusebda i 2 n yibenk n ukayad.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Alemmes n yibenk

Tafelwit s talɣa V d tarrayt n ulemmes n yibenk.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Akersi yettuseggamen s trisiti

Akersi PRIMA s useggem elektrik n teɣzi.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17


Suddes sin n yibenkan

Tawila tunsibt Easyvisus Two tesdukkel sin n yibenkan deg umkan ameẓyan.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Heggi adeg

Azenziɣ s yikaten yettɛawan aheggi n usbeddi d yimussuten n tfelwit.



Easyvisus Two · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Fren talɣa n taggara

Agraw n yiniten n umeskar yeskan ifranen yellan i usdukkel n umkan deg tklinik.



Easyvisus Two · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Ergonomi

Asbeddi yettwaxedmen ilmend n umkan n ukayad.

Ibenkan yettwasuddesen

Ibenkan ttaren amkan-nsen deg yiwet n tayunt tusligt.

Iseɣtayen n umkan

Imussuten yellan seshlen asgunfu n twila.

Tiwuriwin deg usnas

Tayunt n rifraksyun

Iɣbula d yisemliyen

I wakken ad yishil ubrid n umahil

Ṭṭabla iwulmen i yibenkan yettwafrenen

Tamsisɣelt tamagnut 400102 tesɛa tafelwit V n teɣzi ur nettbeddil ara i sin n yibenkan, s umekanizm uslig n uslemmas. Ameskar ur yeskan ara dagi ṭṭabla yettuzzuyun neɣ yettselliden s useggem ibeddi s trisiti. Tajabbuḍt i yettwakkasen, mgal aɣebbar, i uptalmusku d uskiyasku; ibenkan-a ur ttwaḥsaben ara dduklen. Lqem i yizelmaḍen 400102S d afran.

Iɣbula d yisemliyen

Inadiyen d users s yisemliyen

Afus n ṭṭabla yettusekkeṛ s umekanizm deg akk imukan. Anagraw-a yemgarad ɣef usekkeṛ iliktrumagnitik n Curve, RM64 d MS33. Teɣzi s trisiti d tin n ukersi PRIMA amagnu, mačči n tfelwit V. Aslemmas n yibenkan ilaq ad iwulem i wid yettwafrenen deg tesɣunt ASEM. Anadi s uḍar n ukersi yella deg yimallalen n ufran.

Iɣbula d yisemliyen

I useqdec n yal ass

Akersi amagnu d yimallalen n ufran

PRIMA amagnu yesɛa aseggem n teɣzi s trisiti d taduli skay taberkant, war asersi n yiḍarren. Asersi n yiḍarren i yettalayen, tidulin-nniḍen d ubeddel s wanaw-nniḍen n ukersi d ifranen iɣef llan yisemliyen. Leqwamat n ukersi ilaq ad ttwaseɣten deg tesɣunt ASEM. Tawlaft n tazwara tezmer ad tesken asersi n yiḍarren, afus n urifraktyur neɣ imallalen-nniḍen n ufran; askan-nsen ur iḍemmen ara dduklen. Ibenkan d yimallalen ilaq ad ttwasegzin s yinagrawen yemgaraden.

Iɣbula d yisemliyen

Tafawt d trisiti i yibenkan

Atransfurmatur i trisiti n yibenkan yesɛa isemli. Tagejdit n laluminyum tesɛa taftilt LED n tɣuri d tfelwit n users i usenfal n yikayaden. Tabrušurt ur d-tettak ara dagi tazwert neɣ asnagar n trisiti i tayunt meṛṛa; azal n 230 V n yibenk n One ur yettwanɣal ara ihi ɣer Two. Ibenkan i yettaṭṭafen trisiti d yisuturen-nsen ilaq ad ttwaseɣten akked ASEM.

Iɣbula d yisemliyen



Easyvisus Two · Tugna

Iɣbula d yisemliyen

Asbeddi d ufran n twila

Azenziɣ d uslig i temsisɣelt 400102, yeskan ikaten yemgaraden n usbeddi. Ur ddukklen ara deg tsekta timatutin, teɣzi i yettuseggamen neɣ amussu s umutur. Taẓayt ur d-tettunefk ara deg yisebtar-a. Afran n yefus neɣ zelmaḍ, adeg ilelli, akersi d yimallalen ilaq ad ttwasegzin i tklinik. Afus n urifraktyur neɣ asersi n yiḍarren i yettwawalen deg tawlaft ur yettwaḥseb ara d ayen yedduklen s wudem amagnu.

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
F.I.S.O.

Anaw
Easyvisus Two

Tatiknuluji
Tayunt n rifraksyun

Aseqdec akliniki
Itiknulujiwin tusligin

Amḍan n yibenkan
2

Iɣbula d yisemliyen

Teɣzi n ukersi
Aseggem s trisiti; akersi Prima amagnu

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Teɣzi n ṭṭabla
Teɣzi ur nettbeddil ara

Iɣbula d yisemliyen

Asekkeṛ n ṭṭabla
Afus n rifraktyur: asekkeṛ amikanik

Iɣbula d yisemliyen


Amussu n ṭṭabla
Ṭṭabla s talɣa V; amekanizm n uslemmas

Iɣbula d yisemliyen

Tafawt n umkan n umahil
LED n tɣuri


Iɣbula d yisemliyen

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.



Easyvisus Two · Tugna

Iɣbula d yisemliyen



Easyvisus Two · Tugna

Iɣbula d yisemliyen

Asider d yiɣbula

Isemliyen n umeskar

PDF · 15,672 kB

Iɣbula d yisemliyen

- https://www.fiso.info/en/products/refraction_units/easyvisus-two_454/
  Aɣbalu yettwasentem ass n 2026-10-05
- https://www.fiso.info/en/products/refraction_units/easyvisus-two_454/
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-06
- https://www.fiso.info/en/products/refraction_units/easyvisus-two_454/
  Aɣbalu yettwasentem ass n 2026-10-17
- https://www.fiso.info/en/products/refraction_units/easyvisus-two_454/
  Aɣbalu yettwasentem ass n 2026-08-30
- https://www.fiso.info/dynatc/brochure-0497-2d40/frastema_cat-fiso_16luglio_2025_med.pdf
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef Easyvisus Two

D acu Amḍan n yibenkan iɣef llan yisemliyen n Easyvisus Two?

2

Iɣbula d yisemliyen

D acu Teɣzi n ukersi iɣef llan yisemliyen n Easyvisus Two?

Aseggem s trisiti; akersi Prima amagnu — Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anta tawila i ilaq ad sentmeɣ i Easyvisus Two?

Asersi n yiḍarren d yikersiyen-nniḍen n ufran ur ttwaḥsaben ara ttwasebdedden.

Iɣbula d yisemliyen

Anda ara afeɣ isemliyen n Easyvisus Two?

Isemliyen n umeskar (PDF) — https://asem-medical.com/storage/catalogue-enrichment/fiso-20261006/fiso-original-2025.pdf

Iɣbula d yisemliyen

Anwa talɣa d teɣzi n tfelwit i d-yettwasegzan?

Easyvisus Two 400102 yessemras tafelwit s talɣa V i sin n yibenkan, s teɣzi tuqfilt d umekanizm n useggem n yibenkan deg tlemmast.

Iɣbula d yisemliyen

Anwa anaw n usekkeṛ i d-yettwasegzan?

Tawriqt teskan asekkeṛ amikanik n uɣil n tfelwit deg yal adeg. Asekkeṛ aliktrumagnitik n tayunin-nniḍen ur as-yettunefk ara.

Iɣbula d yisemliyen

Asbeddi azelmaḍ yettwahegga?

Tawriqt teskan lqem azelmaḍ afrayan 400102S. Akersi d yimerna ilaq ad ttusegzin deg ufran n twila.

Iɣbula d yisemliyen

Ifarisen icudden



Panvisus Rondo

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt Panvisus Rondo i sin n yibenkan, ilmend n twila.

Wali afaris



Easyvisus Curve

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 2 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 2 n yibenk n ukayad.

Wali afaris



Easyvisus MS33 — 3 n yibenk

F.I.S.O.

Tayunin n usekyed n tɣuri

Tayunt n rifraksyun i 3 n yibenk n ukayad.

Wali afaris

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Tetiknulujiyin tusliḍin

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: ad706d73d80f