Tiftilin s tneḍrit

Elara 900

Taftilt s tecraḍ, s umikruskuḍ agalili d 5 n yisemeɣriyen.

Taɣara Haag-Streit. S usdukkel n tiḍwant Haag-Streit d tergonomi tanernant, Elara 900 tettak abrid amaynut i ukayad n wallen. Tamuɣli taklinikant: asenqed n yibeddilen imeẓyanen n tiṭ s usezri n tafat ɛlayen d unagraw aḍwan yedduklen, yettwasnasen, iṣeḥḥan. Lehna yemden: ergonomi i umassan d useqdec yelhan i umuḍin, s uɛeyyu n twali imeẓyanen d yikayaden iruraden. Abrid yeshilen: inadiyen isehlen, iseɣtayen imiren kan, i urbaḥ n wakud n yal ass. Taflest taswisrit: lqewwa, taɣara yettwasenqden d uṣebber i useqdec s waṭas. Asidref n zdat: ibenk i yezemren ad yennerni, i yeḍfaren itiknulujiwin d yiberdan imaynuten n tdawsa n wallen. Elara 900: telqay, lehna d wakud ɣezzifen i yikayaden n tura d zdat.



Elara 900 · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-10

Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Tafat s usnay n tecraḍ

Tafat s usnay i tmugli n tɣessiwin n tiṭ.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Asenqed s uḍaw infraruz

Aḍaw infraruz i usenqed n tɣuddanin Meibomius.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Asenqed s umutur

Asenqed s umutur n teɣzi n yibenk d usemɣer.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Tiḍwant Haag-Streit

Anagraw aḍwan yettwaxedmen i twali s telqay n uḥric n zdat.

Asers yettwaḥkamen

Iseggimen ttaɛawanen asers yeshilen n yal ass.

Imussuten isehlen

Inadiyen qqimen ttwaseggmen ɣef ukayad akliniki.

Tiwuriwin deg usnas

Wali s telqay

Anagraw aḍwan yettɛawan akayad s telqay n lqerniya d uḥric n zdat.

Iɣbula d yisemliyen

Af-d akna iwulmen

Taɣessa n umkan teshel iseɣtayen ilaqen i twali.



Elara 900 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Ḥkem tafawt

Anagraw n tafawt yettak twennaḍt iwulmen i useqdec akliniki.



Elara 900 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Seggem tafat ilmend n twali

Aɣbalu n usufeɣ P-Type yettak aṭas n yiniten d tezɣiwin n tafat. Amsizdeg awraɣ s umutur yettnerni ɣer yiseggamen n twali.



Elara 900 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Semmus n yiswiren n usemɣer

Amikruskuḍ agalili yeɣḍel 6.3×,10×,16×,25× d 40×. Asemɣer yettuseggam seg joystick s unadi s umutur.



Elara 900 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Taftilt n tfexsiḍt i d-yettwasnulfan tikkelt-nniḍen

Tamuɣli taklinikant s telqay — senqed yal aferdis s taflest

Taflest taswisrit

Yezmer ad yennerni — yettwaxdem i tura d zdat

Xdem s tnezzut

Ikar iwurmanen d tergonomi iwulmen i lehna n umassan d umuḍin.

I wakken ad yishil ubrid n umahil

Ales afran n ubrid n ukayad

Iseggamen i yellan si tazwara sdukklen tafat, amsizdeg, tafexsiḍt d usemɣer. Tazelqact tettak-d abrid n ukayad war ma nɛawed yal aseggem s ufus.



Elara 900 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Ṭṭef sakin sismel tugniwin d tvidyutin

Ibenk yettaṭṭaf tugniwin d tvidyutin seg ugrudem-is. Asiweḍ ɣer EyeSuite neɣ akaram n umuḍin DICOM yeqqen ɣer tawennaḍt d turagin yettwafernen.

Tuqqna DICOM teḥwaǧ turagt; senqed asdukkel i d-yettunefken.

Iɣbula d yisemliyen

I useqdec n yal ass

Seggem adeg seg umkan n umahil

Teɣzi s umutur, lentes n twali yenḍan d tmessendet n lḥenk yettuseggamen ttɛawanent aseggem n wadig. Inadiyen igejdan duklen ɣer joystick.

Iɣbula d yisemliyen



Elara 900 · Tugna

Iɣbula d yisemliyen

Fren tugniwin d tonometri

Anṭaf astiriyuskupik 3D yeḥwaǧ turagt-is. Tonometri Goldmann teḥwaǧ AT 900 n wanaw R imṣada, i d-yettrennun ɣer teftilt.

Ifran d wayen yellan deg yal tamnaḍt ilaq ad ttwasentmen; ulac amerna i nettwaḥseb d ayen yeddan.

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
Haag-Streit

Anaw
Elara 900

Tatiknuluji
Taftilt n tfexsiḍt

Aseqdec akliniki
Adiyagnustik d biyumetri

Asimɣuren
6.3× · 10× · 16× · 25× · 40×

Iɣbula d yisemliyen

Amḍan n yisimɣuren
5

Iɣbula d yisemliyen

Tiḍwant n umikrusku
Agali

Iɣbula d yisemliyen

Tisekta
560 × 401 × 330 mm

Iɣbula d yisemliyen

Taẓayt
15.2 kg

Iɣbula d yisemliyen



Elara 900 · Tugna

Iɣbula d yisemliyen



Elara 900 · Tugna

Iɣbula d yisemliyen

Asider d yiɣbula

- https://haag-streit.com/en/products/categories/general-diagnostics/slit-lamps/slit-lamps
  Aɣbalu yettwasentem ass n 2026-10-04
- https://us.haag-streit.com/Campaigns/elara/900
  Aɣbalu yettwasentem ass n 2026-10-17
- https://us.haag-streit.com/Campaigns/elara/900
  Aɣbalu yettwasentem ass n 2026-08-30
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/a-superior-clinical-view-uniq-68bfd159b9f2e
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/a-superior-clinical-view-uniq-68bfd159b9f2e
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/enhanced-efficiency-uniq-68b879601e8ae
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/enhanced-efficiency-uniq-68b879601e8ae
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/optimized-ergonomics-uniq-68b8798975362
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/faq
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/faq
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/enhanced-efficiency-uniq-68b879601e8ae
  Aɣbalu yettwasentem ass n 2026-10-18
- https://us.haag-streit.com/products/categories/general-diagnostics/slit-lamps/slit-lamp-elara-900/faq
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef Elara 900

D acu Asimɣuren iɣef llan yisemliyen n Elara 900?

6.3× · 10× · 16× · 25× · 40×

Iɣbula d yisemliyen

D acu Amḍan n yisimɣuren iɣef llan yisemliyen n Elara 900?

5

Iɣbula d yisemliyen

Tugniwin 3D ddant-d s tazwara?

Ala. Ameskar yenna belli turagt 3D ilaq; tawila i d-yettunefken ilaq ad tt-tessegzi.

Iɣbula d yisemliyen

Tuqqna DICOM teḥwaǧ turagt?

Ih. Asebter n umeskar yeqqen asdukkel n DICOM worklist ɣer turagt Connectivity icudden.

Iɣbula d yisemliyen

Anwa tonometr Goldmann i d-yettwabedren i Elara 900?

Asebter n umeskar yeskan AT 900 n wanaw R; ilaq ad yettwasegza ma yedda-d deg usiweḍ.

Iɣbula d yisemliyen

Ifarisen icudden




BI 900

Haag-Streit

Tiftilin s tneḍrit

Taftilt s tecraḍ, amikruskuḍ Kepler, 2 n yisemeɣriyen d LED.

Wali afaris



BM 900

Haag-Streit

Tiftilin s tneḍrit

Taftilt s tecraḍ, amikruskuḍ Kepler d yisemeɣriyen 10×/16×.

Wali afaris



BQ 900

Haag-Streit

Tiftilin s tneḍrit

Taftilt s tecraḍ, amikruskuḍ agalili, 5 n yisemeɣriyen d LED.

Wali afaris



HS-5000

Huvitz

Tiftilin s tneḍrit

Taftilt s tecraḍ d tafat n tburt, s tḍwant tagalilit.

Wali afaris

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Asenqed d biyumetri

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: 0c773c175db6