Aktazal n wurti n yiẓri

Octopus 600

Aperimetr s perimetri tusbiḍt SAP d perimetri Pulsar.

Tulmisin titiknikanin tigejdan 1. Ikayaden imagnuten n pirimetri ur nettḥerrik ara: • Ikayaden n umnar amellal ɣef umellal i uḥric alemmas n twali. • Tarrayt TOP taruradt: seg 2 ar 4 n tesdatin i ukayad yemden n umnar. • Tagrumma n yikayaden i yisuturen iklinikiyen yemgaraden: Aḥric alemmas: G, 32, 30-2, 24-2. Ikayaden n makula: M (yebnan ɣef tfizyuluji) d 10-2. • Tuqqna akked talɣiwin n usiggez Octopus d HFA. 2. Asenqed arurad: • TPS (Tendency Perimetry Screening): Yessemgarad gar uḥric n twali amagnu d win ur nelli d amagnu deg ddaw n yiwet n tesdat. Isfuklen usligen ttwaskanen 3 n tikkal i usentem n yigmaḍ. • Pirimetri Pulsar: Asenqed uzwir n glukuma. 3. Tasleḍt s telqay n yigmaḍ: • Tesdukkel lixṣas n twali ilmend n yifilan inerviyen n retin. • Tiḥulfa ɛlayen i ubeddel azwir n glukuma. • Aseɣẓan EyeSuite — Tasleḍt n usfari: Yeskan abeddel adigan neɣ yezzuzren d urured-ines. Iskanen izegzayen s yisfuklen d yizamulen isehlen. 4. Ergonomi d tuqqna n yiseɣẓanen: • Aseɣẓan EyeSuite: Yesdukkel igmaḍ n waṭas n yibenkan Haag-Streit. Tuqqna akked yinagrawen EMR war aseɣẓan wis kraḍ. • Ibenk ameẓyan d ufessas i yimukan imeẓyanen. • Yettwaseqdec deg tklinikin s waṭas n yimuḍan. 5. Ifranen n uskan n yisefka: • Asiggez amagnu 7-deg-1: igmaḍ s talɣiwin yemgaraden i zemren ad ttuseggment. • Askan deg uẓeṭṭa: anekcum seg PC s tuqqna EyeSuite.



Octopus 600 · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

> Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Perimetri tusbiḍt n tlemmast

Ikayaden n wurti n tlemmast, amellal ɣef umellal.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Tasleḍt tamatut n unerni

Asenqed n yikayaden imseḍfaren d unerni-nsen amatu.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17


War tɣerɣert taberkant

Ikayaden zemren ad ilin berra n tɣerɣert taberkant.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

Wali aḥric alemmas

Iwahilen usligen ssudusen akayad n uḥric alemmas n twali.



Octopus 600 · Tugna


Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Asenqed s tnezzut

Abrid n usenqed yettnhu anṭaf, yettak agmuḍ i yettwaseqdacen imiren kan.



Octopus 600 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Isemliyen n uḍfar

Tuqqna EyeSuite tesuddes igmaḍ i tɣuri d userwes.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Aḥric n twali

Iwahilen usligen i ukayad n twuri n tiṭ.

Abrid yettwanhan

Anermis yesuddes timellilin n ukayad.

Tasleḍt n uḍfar

Igmaḍ ttaɛawanen aserwes deffir wakud.

Tiwuriwin deg usnas

Pirimetri SAP ur nettḥerrik ara

Yemṣada i Basic / Pro; TOP SAP yeqqen ɣer lqem.

Iɣbula d yisemliyen


Pirimetri Pulsar

Yemṣada i Basic / Pro; TOP SAP yeqqen ɣer lqem.

Iɣbula d yisemliyen

Ktil aḥric alemmas n twali

Octopus 600 ixeddem perimetri tabeddant amellal-ɣef-umellal n uḥric alemmas. Ihilen G, M, 24-2 d 10-2 qqnen ɣer yibḍanen yemgaraden n yireqqiḍen.



Octopus 600 · Tugna


Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Frez gar usenqed amezwaru d wemnar

Asenqed arurad d win n tɣara. Tarrayt Pulsar tessemras aseḥrek n twali s umgired n tafat d tifrat yettbeddilen; ur tettarra ara deg umkan n perimetri tameslugant n uḍfar.



Octopus 600 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Asenqed n uṭṭef n tmuɣli

Isenqed s wudem awurman tuččiwin n uṭṭef n tmuɣli.

Tuččiwin n uṭṭef n tmuɣli ssebba-nsent d ur nṣebbir ara n umuḍin, zemrent ad d-fkent aḥric n twali ur neṣḥiḥ ara. Asenqed Haag-Streit yettḥebbis akayad ma tmuɣli ur tessewseɛ ara, syin yettkemmil s wudem awurman. Aya yettɛawan i wakken yal agaz n uḥric n twali ad yettwakayad.

Aseqdec deg twennaḍt yecɛelen

• Agdil TFT yeskan isfuklen ur nettbeddil ara s tafawt n twennaḍt.
• Yettwaseqdec berra n texxamt taberkant, war isuturen usligen n tafat.

Ergonomi d useqdec yeshilen

• Taḥnunt n tririt iwulmen:
Asentem s ṣṣut d umussu yettwasnerna.
Taɣessa yeshilen i tririt deg yimukan yemgaraden.
• Abrid n umahil yettwanhan:
Yettɛawan amsekyad ad iḍfer timellilin ilaqen i yikayaden iṣeḥḥan.
Yeshel almad d usebɣer n tɣara n yigmaḍ.

Asers yeshilen n umuḍin

• Tilentiyin n ukayad timeqqranin:
Seshilent asers s umekanizm amagnitik.
• Agdal n tiṭ n tafat ur neččur ara:
Yettḥebbis isfuklen i tiṭ iɣelqen, maca yeğğa-yas tafat tafessast.
Ur yeḥwağ ara taduli n tiṭ; yesrus ayen ur nettwagdalen ara s lekmal, ur yeḥwağ ara asezdeg s ddwa.
• Asers yelhan:
Asers aknun kra ɣer zdat d ddaw yeshil ugar.
Asersi n tmart yella d afran i ugar n lehna.

Isnasen

• Iwulem i tklinikin n wallen i yeḥwağen ibenk ameẓyan, yettwasnasen, yeshilen i useqdec.
• Yettwaxdem i yimukan anda asehel n useqdec d tɣara n yigmaḍ d ayen agejdan.

Iferdisen igejdan

Asenqed arurad:

Yessemgarad gar uḥric n twali amagnu d win ur nelli d amagnu deg ddaw n yiwet n tesdat s pirimetri Pulsar.

Tasleḍt s telqay n glukuma:

Askan n yikayaden n uḥric n twali s yidisan yemgaraden, s tesleḍt n yigrawen d tesleḍt apular i userwes n tɣessa.

Ergonomi d useqdec yeshilen:
Ibenk ameẓyan i yettwaseqdacen berra n texxamt taberkant, s users iwulmen i umsekyad d umuḍin.

Tuqqna n yiseɣẓanen:

EyeSuite yeshel tuqqna ɣer ubrid n umahil n tklinik, s tuqqna EMR d usuddes n yisefka deg yiwen n wadeg.

I wakken ad yishil ubrid n umahil

Senqed aṭṭaf n twali d wadig


Blink Control d Pupil Position Control senqaden aṭerref d tuffɣa n tlemmast. Dart Control d Contact Control ttemmaden asenqed n tiṭ d wadig n umuḍin.

Iɣbula d yisemliyen

Ɣer azrar n yiḥricen n twali

EyeSuite yeskan iḥricen iman-nsen d ubeddel. Tisleḍin Cluster Trend d Polar Trend d timesluganin deg Pro, d ifran deg Basic ilmend n tfelwit tamirant.

Tiwuriwin ilmend n Basic/Pro d yifran yermden.



Octopus 600 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

I useqdec n yal ass

Xdem berra texxamt taberkant

Agdil TFT yettağğa aseqdec berra texxamt taberkant deg waddaden n umeskar. Agdil amḍufan, anadi n tririt d lentes n usekyed ttɛawanen akayad.

Iɣbula d yisemliyen

Fren tiwuriwin n umkan

TOP, tisleḍin n ubeddel, DICOM d unekcum ɣer taffa yebɛden ilaq ad ttwasenqden i Basic neɣ Pro. Ulac azegrir n perimetri s umussu i d-yettwabedren i 600.

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
Haag-Streit

Anaw
Octopus 600

Tatiknuluji
Apirimetr

Aseqdec akliniki
Adiyagnustik d biyumetri

Tisekta n yisfuklen
SAP Goldmann III · Pulsar 5°

Iɣbula d yisemliyen

Aḥric afellay n ukayad
30 °

Iɣbula d yisemliyen

Tuqqna
Ethernet · EMR

DICOM d afran, ur yettwaḥseb ara yeddukel.

Iɣbula d yisemliyen


Tisekta
467 × 508 × 500 mm

Iɣbula d yisemliyen

Taẓayt
12.7 kg

Iɣbula d yisemliyen


DICOM d afran, ur yettwaḥseb ara yeddukel.



Octopus 600 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen



Octopus 600 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Asider d yiɣbula

- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20600/Brochures/English/HS_octopus_600_br_7220354_perimetry_simplified.pdf
  Aɣbalu yettwasentem ass n 2026-10-05
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-600
  Aɣbalu yettwasentem ass n 2026-10-17
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-600
  Aɣbalu yettwasentem ass n 2026-08-30
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20600/Brochures/English/HS_octopus_600_br_7220354_perimetry_simplified.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20600/Brochures/English/HS_octopus_600_br_7220354_perimetry_simplified.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20600/Brochures/English/HS_octopus_600_br_7220354_perimetry_simplified.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-600
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20600/Brochures/English/HS_octopus_600_br_7220354_perimetry_simplified.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-600
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-600
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20600/Brochures/English/HS_octopus_600_br_7220354_perimetry_simplified.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-600
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef Octopus 600

D acu Tisekta n yisfuklen iɣef llan yisemliyen n Octopus 600?

SAP Goldmann III · Pulsar 5°

Iɣbula d yisemliyen

D acu Aḥric afellay n ukayad iɣef llan yisemliyen n Octopus 600?

30 °

Iɣbula d yisemliyen

Anta tawila i ilaq ad sentmeɣ i Octopus 600?

DICOM d afran, ur yettwaḥseb ara yeddukel.

Iɣbula d yisemliyen

Octopus 600 ixeddem perimetri s umussu?

Ala. Tafelwit Basic/Pro ur d-tebdir ara askar s umussu i wanaw-a.

Iɣbula d yisemliyen

Asenqed arurad yettawi amkan n ukayad ɣef wemnar?

Ala. Asenqed i yettwasegzan yerza taɣara; ikayaden n wemnar d usefhem-nsen qqimen yemgaraden.

Iɣbula d yisemliyen

Tarrayt TOP tedda-d deg yal lqem?

Tafelwit tamirant teskan-itt d afran deg Basic, d tameslugant deg Pro.

Iɣbula d yisemliyen

Ifarisen icudden



Octopus 900

Haag-Streit

Aktazal n wurti n yiẓri

Perimetri ur nettbeddil ara; tiwuriwin-nniḍen ilmend n twila.

Wali afaris



HVF-100

Huvitz

Aktazal n wurti n yiẓri

Akat n wurti n tmugli s perimetri ur nettbeddil ara, amellal ɣef umellal.

Wali afaris

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Asenqed d biyumetri

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: c0b86784bf6b