Aktazal n wurti n yiẓri

# Octopus 900

## Perimetri ur nettbeddil ara; tiwuriwin-nniḍen ilmend n twila.

Octopus 900 d apirimetr s tqubbet n Haag-Streit i ukayad n uḥric alemmas d yiri n twali. Tinasant tesɛa pirimetri ur nettḥerrik ara d yikayaden imussuten ilmend n twuriwin iremden. EyeSuite yesuddes igmaḍ, ikayaden izwuran d tesleḍt n usfari. Basic d Pro yemgaraden s yigrawen d yituragen i d-yettunefken; TOP, tarrayin tusligin, amussu Goldmann d yinermisen ilaq ad ttwaseɣten i twila ASEM, war leḥsab akk tzemmar n tewsit dduklen.



Octopus 900 · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-09

> Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

# Iferdisen igejdan

## Perimetri n wurti ummid

Ikayaden n wurti n tlemmast d yiri s yinagrawen usbiḍen.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

## Akayad n wurti s umussu

Ikayaden s umussu ilmend n lqem neɣ turagt.

Tafrant ilmend n lqem

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17


## Tasleḍt tamatut n unerni

Tasleḍt n unerni s yikayaden imseḍfaren n wurti n tmugli.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-17

## Kayad akk aḥric n twali

Iwahilen iɣef llan yisemliyen sdukkulen ikayaden n tlemmast d yiri deg yiwen n unagraw.




Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

### Seggem tasleḍt ɣer yiwet n tarrayt

EyeSuite yesdukkel isemmal d yiskanen iwulmen i ubrid n pirimetri.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

### Ḍfer abeddel

Tugniwin n usfari seshlent aserwes n yikayaden n twuri deffir wakud.



Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

### Aḥric n twali

Iwahilen usligen i ukayad n twuri n tiṭ.

### Abrid yettwanhan

Anermis yesuddes timellilin n ukayad.

### Tasleḍt n uḍfar

Igmaḍ ttaɛawanen aserwes deffir wakud.

## Tiwuriwin deg usnas

### Pirimetri ur nettḥerrik ara

Amussu / TOP / SWAP qqnen ɣer Basic / Pro d yifrangen; ur ttwaḥsaben ara dduklen.

Iɣbula d yisemliyen

## Wali aḥric alemmas d yiri n twali

Aḥric n yiri: 180°. Amasay n tqubbet: 30 cm. Tisekta n yisfuklen Goldmann: I / II / III / IV / V. Tilisa timussanin: 47 dB. Asfukel yettwasken deg tqubbet. Iwahilen iɣef llan yisemliyen sdukkulen 24-2 d 10-2 i tlemmast, G d M d yikayaden n yiri. Amassan yefren ahil ilmend n usteqsi akliniki; ulac afran awurman n udiɣnustik neɣ n tzemmar tusrifin n unhar i d-yettwakkasen seg tilin n ukayad kan.

Iɣbula d yisemliyen


## Tarrayin d yibriden: senqed Basic neɣ Pro

Pirimetri amagnu amellal/amellal yesɛa isemliyen deg Basic d Pro. Tafelwit tessemgarad igrawen d ifranen n Basic i TOP, Goldmann imussu, anili/awraɣ d kra n yikayaden usligen. Akud n tarrayin i d-yettunefken d TOP 2–4 n tesdatin, Dynamic 6–8 n tesdatin d Normal 10–12 n tesdatin; maċċi d leḍman n wakud n yal akayad. Amussu d yifrangen iremden s tidet ilaq ad ttwasegzin deg tesɣunt.

Iɣbula d yisemliyen

## Seg tlemmast ɣer yiran

Taqubbat Goldmann teɣḍel ikayaden n tlemmast d yiran. Tawriqt teskan aḥric ammuɣ n 180° d yihilen gar-asen G, M d Esterman n snat n wallen.




Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

## Suneɣ izuptir s umussu

Perimetri Goldmann s umussu tettak afran n yivikturen, izuptir d yileqman n ukayad. Armad-is yeqqen ɣer lqem d yifran yettwasbedden.

D afran deg Basic; d ameslugan deg Pro ilmend n tfelwit tamirant n umeskar.

Iɣbula d yisemliyen

# I wakken ad yishil ubrid n umahil

## Asenqed n uṭṭef n tmuɣli d tɣuri n yigmaḍ

Tawsit teskan Blink Control, Pupil Position Control d usfuɣel awurman AET. EyeSuite yesdukkel tasleḍt tamatut d tugniwin n usfari; Cluster / Polar Trend qqnen ɣer yigrawen n twila. Asenqed n tmuɣli d yisemmal ttaɛawanen asefhem n umassan war abeddel n usteqsi-ines ɣef taflest n yal akayad.

Iɣbula d yisemliyen

## Tuqqna ɣer uẓeṭṭa d yidufa n yimuḍan

Ethernet d twennaḍt EyeSuite ssudusen ikayaden d yimukan n ukayad. Tabrušurt tebder akter seg Octopus iqburen d yisefka HFA iwulmen. DICOM d afran; inermisen EMR, leqwamat d yisquḍḍen ilaq ad ttwaseɣten deg usbeddi s tidet. Ulac ttuɛud n tuqqna akked akk idufa n yimuḍan.

Iɣbula d yisemliyen


## Ɛiwen abrid n ukayad

Abrid yettusegzan, tinadin i umuḍin s waṭas n tutlayin d usenqed n uṭṭaf n twali sseggamen akayad. Iseḥriken n twali i iruḥen seld aṭerref zemren ad d-ttwasknen daɣen.

Iɣbula d yisemliyen

## Serwes iḥricen n twali s wakud

EyeSuite yesdukkel iḥricen n twali, tisleḍin timatutin d tmuɣliwin n ubeddel. Iskanen Octopus neɣ HFA ttɛawanen taɣuri ilmend n tannumi n ukabini.



Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

# I useqdec n yal ass

## Senqed aseɣti n usekyed

Aseḍfar n lente s umutur d uḍfar n umeccaq-is ddant deg twila yettwasegzan. Asenqed n tberract d uṭṭaf n twali ttɛawanen aseggem n wadig.



Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

## Asbeddi d twila

Tehri: 648 mm. Telqay: 519 mm. Teɣzi: 796 mm. Taẓayt: 25 kg. Asnagar: 50/60 Hz. Unaḍ yeskan snat n tilisa n trisiti: 100–120 V AC neɣ 220–240 V AC. Aseqdec yettwammelen d 145 / 165 VA ilmend n twila. Amkan, asersi, ituragen d yibenkan n uẓeṭṭa ilaq ad ttwaseɣten akked ASEM. Taẓayt n usiweḍ n 40 kg s ubukel maççi d taẓayt n yibenk; ur tettwaseqdec ara d taẓayt n yibenk iman-is.

Iɣbula d yisemliyen

# Fren ayen yellan deg Basic neɣ Pro

SWAP, Flicker, iseḥriken izeggaɣen, askar s umussu d TOP ur ddin ara deg yal lqem. Asuter n ssuma ilaq ad yessegzi lqem, ikayaden, tisleḍin d yigrudmawen yellan.

Iɣbula d yisemliyen

# Isallen itiknikanen d yisemliyen

### Ticreḍt
Haag-Streit

### Anaw
Octopus 900

### Tatiknuluji
Apirimetr

### Aseqdec akliniki
Adiyagnustik d biyumetri

### Tisekta n yisfuklen
Goldmann I · Goldmann II · Goldmann III · Goldmann IV · Goldmann V

Iɣbula d yisemliyen

### Tazmert tafellayt n usfukel
10000 asb


Iɣbula d yisemliyen

### Aseɣẓan
EyeSuite

Iɣbula d yisemliyen

### Taẓayt
12.7 kg

Iɣbula d yisemliyen




Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen



Octopus 900 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

# Asider d yiɣbula

## Isemliyen n umeskar

PDF · 10,177 kB

Iɣbula d yisemliyen

## Isemliyen n umeskar

PDF · 5,010 kB

Iɣbula d yisemliyen

- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-05
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-06
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20900/Brochures/English/HS_Oktopus_900_br_7220444_english.pdf
  Aɣbalu yettwasentem ass n 2026-10-06
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20900/Instructions%20for%20use/English/1500_7220343_04121_IFU_OCTOPUS_900_01_en_web.pdf
  Aɣbalu yettwasentem ass n 2026-10-06
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-17
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-08-30

- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20900/Brochures/English/HS_Oktopus_900_br_7220444_english.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20900/Brochures/English/HS_Oktopus_900_br_7220444_english.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20900/Brochures/English/HS_Oktopus_900_br_7220444_english.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/2%20Products/Speciality%20diagnostics/Perimetry/Octopus%20900/Brochures/English/HS_Oktopus_900_br_7220444_english.pdf
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-18
- https://haag-streit.com/en/products/categories/specialty-diagnostics/perimetry/octopus-900
  Aɣbalu yettwasentem ass n 2026-10-18

# Isteqsiyen-nwen ɣef Octopus 900

### D acu Tisekta n yisfuklen iɣef llan yisemliyen n Octopus 900?

Goldmann I · Goldmann II · Goldmann III · Goldmann IV · Goldmann V

Iɣbula d yisemliyen

### D acu Tazmert tafellayt n usfukel iɣef llan yisemliyen n Octopus 900?

10000 asb

Iɣbula d yisemliyen

## Anda ara afeɣ isemliyen n Octopus 900?

Isemliyen n umeskar (PDF) — https://asem-medical.com/storage/catalogue-enrichment/haag-streit-remaining-20261006/octopus900-original.pdf Isemliyen n umeskar (PDF) — https://asem-medical.com/storage/catalogue-enrichment/haag-streit-remaining-20261006/octopus900-ifu-original.pdf

Iɣbula d yisemliyen

## Askar Goldmann s umussu yedda-d yal tikelt?

Ala. Tafelwit tamirant teskan-it d afran Basic, d tawuri tameslugant Pro.

Iɣbula d yisemliyen

## Nezmer ad nseqdec askan n wanaw HFA?

Ih. EyeSuite yettağğa afran gar iskanen Octopus d HFA i yiḥricen n twali d yineqqisen i d-yettwasegzan.

Iɣbula d yisemliyen

## Akayad Esterman n snat n wallen yettwasegza?

Ih. Yella gar yihilen i d-yettwasknen i uktil n uḥric n twali n snat n wallen.

Iɣbula d yisemliyen

# Ifarisen icudden




## Octopus 600

Haag-Streit

Aktazal n wurti n yiẓri

Aperimetr s perimetri tusbiḍt SAP d perimetri Pulsar.

Wali afaris



## HVF-100

Huvitz

Aktazal n wurti n yiẓri

Akat n wurti n tmugli s perimetri ur nettbeddil ara, amellal ɣef umellal.

Wali afaris

# Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Asenqed d biyumetri

# Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

### Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

### Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

### Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

### Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: 39abe07ac1db