Asenqed amezwaru n yiẓri

HVS-1

Asenqed n rifraksyun n snat n wallen seg lbeɛd.

Asenqed arurad n twali i yilufanen d warrac imeẓyanen. Akat deg ddaw n yiwet n tasint. Ikaten iṣeḥḥan ugar s tnezmart n LED yettwasnernan, d useqdec yeshilen i umseqdac. Mmeslay.



HVS-1 · Tugna

Iɣbula d yisemliyen

Suter tasɣunt n ssuma

Wali afaris

Tettwahegga ass n 2026-11-09

> Taferkit-a n ufaris thegga-tt ASEM. Ur tettbeddil ara isemliyen unsiben n umeskar neɣ tinaḍin n useqdec.

Iferdisen igejdan

Asenqed n snat n wallen

Akat n rifraksyun n snat n wallen seg lbeɛd.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-05

Azilal n rifraksyun

SPH seg −8 ar +8 D; CYL alamma 3 D.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-10-05

Asenqed yeshilen

Asers n snat n wallen d ubrid yettwanhan ssudusen asenqed n twali.



HVS-1 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Ɣer igmaḍ imiren kan

Agdil yesdukkel ikaten igejdan d ussefhem-nsen i usenqed.

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Azen ɣer HIIS

Isefka iɣef llan yisemliyen zemren ad ddun deg ubrid Huvitz i uḍfar deg yiwen n wadeg.



HVS-1 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen · Aɣbalu yettwasentem ass n 2026-08-30

Akat awurman

Abrid uslig i unṭaf n yisefka n rifraksyun.

Isefka yettwasuddesen

Isallen iwulmen qqimen yellan deg unermis.

Abrid yeqqnen

Igmaḍ zemren ad ddukklen deg ubrid n rifraksyun n tklinik.

Tiwuriwin deg usnas

Asenqed n rifraksyun n snat n wallen seg lebɛid

Iɣbula d yisemliyen

Kayad tamuɣli seg waggug war asersi n tamart

HVS-1 ixeddem akayad amezwaru n snat n wallen seg wazal n 1 m, s wakud n uktil ara yilin ar yiwet n tasint. Talɣa-ines i yettwattawin tettağğa akayad n uqcic i yeṭṭef win it-iɛawnen.

- Aktil n snat n wallen seg wazal n 1 m.
- Akud n uktil yettwabder: ar yiwet n tasint.



HVS-1 · Tugna

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Arifraktumitr yettwasiwaḍen i warrac

Telqay ɛlayen d useqdec deg yimukan yemgaraden, ladɣa ikayaden deg umussu neɣ tdawsa n warrac, s unermis yeshilen.

Rifraksyun

Asfir.

• Tilisa: -8.0D ar +8.0D.
• Isurifen: 0.01D, 0.25D.
• Telqay:
Seg -3.50D ar +3.50D: ±0.5D.
Seg -8.0D ar < -3.50D: ±1.0D.
Seg +3.50D ar +8.0D: ±1.0D.

Asilindr.

• Tilisa: -3.0D ar +3.0D.
• Isurifen: 0.01D, 0.25D.
• Telqay:
Seg -1.50D ar +1.50D: ±0.5D.
Seg -3.0D ar < -1.50D: ±1.0D.
Seg +1.50D ar +3.0D: ±1.0D.

Agellus.

• Tilisa: 0° ar 180°.
• Isurifen: 1°.
• Telqay: ±10°.
• Tawtilt: i wazalen n usilindr (CYL) > 0.50D.

Ikaten n tpupilt n tiṭ

Tasekta n tpupilt n tiṭ.

• Tilisa: 4.0 mm ar 9.0 mm.
• Isurifen: 0.1 mm.
• Telqay: ±0.4 mm.

Amecwar gar tpupilin n wallen.


• Tilisa: 35 mm ar 80 mm.
• Isurifen: 1 mm.
• Telqay: ±1.5 mm.

I wakken ad yishil ubrid n umahil

Ṭṭef tafat i d-terra retina

Tatiknulujit Focus LED d tesleḍt n tafat i d-terra retina ttaken isefka n ukayad amezwaru n rifraksyun d usuddus n wallen i d-yettbanen deg ugdil.

- Tasleḍt n tafat i d-terra retina.
- Asken n yisefka n rifraksyun d usuddus.



HVS-1 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

I useqdec n yal ass

Sekyed igmaḍ s uẓeṭṭa

Ameskar isbadu amṣada DICOM d usekyed s PC akked HIIS-1. Tiwuriwin-a snernayent asekyed n yisefka n ukayad amezwaru berra n ugdil n yibenk afus.

Ilmend n twila

- Amṣada DICOM i d-yesbadu Huvitz.
- Asekyed s PC deg twennaḍt HIIS-1 yettuseggmen.



HVS-1 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Iɣbula d yisemliyen

Isallen itiknikanen d yisemliyen

Ticreḍt
Huvitz

Anaw
HVS-1

Tatiknuluji
Arifraktyur / arifraktyur awurman

Aseqdec akliniki
Adiyagnustik d biyumetri

Tazwert taddayt n usfir
-8 D

Iɣbula d yisemliyen

Tazwert tafellayt n usfir
8 D

Iɣbula d yisemliyen

Tazwert tafellayt n usilindr
3 D

Iɣbula d yisemliyen

Akdu adday n tpupilt n tiṭ
4 mm

Iɣbula d yisemliyen

Tuqqna
USB 2.0 · Wi-Fi 802.11 b/g/n

Iɣbula d yisemliyen

Aseɣẓan
HIIS-1; tuqqna DICOM tettwammel

Iɣbula d yisemliyen

Asaɣ n trisiti

15 V DC; 2 A

Iɣbula d yisemliyen

Taẓayt

1 kg

Iɣbula d yisemliyen



HVS-1 · Tugna

Iɣbula d yisemliyen



HVS-1 · Tugna

Tugna n umeskar. Tawila yettwaskanen tezmer ad temgarad ɣef tin ara d-yettunefken.

Iɣbula d yisemliyen

Asider d yiɣbula

- https://huvitz.com/wp-content/uploads/2024/09/20250319_HVS-1_kor_catalog_v3.pdf
  Aɣbalu yettwasentem ass n 2026-10-05
- https://huvitz.com/en/product/hvs-1/
  Aɣbalu yettwasentem ass n 2026-10-05
- https://huvitz.com/en/product/hvs-1/
  Aɣbalu yettwasentem ass n 2026-08-30
- https://huvitz.com/en/product/hvs-1/
  Aɣbalu yettwasentem ass n 2026-10-18

Isteqsiyen-nwen ɣef HVS-1

D acu Tazwert taddayt n usfir iɣef llan yisemliyen n HVS-1?

-8 D

Iɣbula d yisemliyen

D acu Tazwert tafellayt n usfir iɣef llan yisemliyen n HVS-1?

8 D

Iɣbula d yisemliyen

Akayad amezwaru yettaɣ amkan n rifraksyun tummidt?

Ameskar iskan HVS-1 d ibenk n ukayad amezwaru n tmuɣli. Igmaḍ-is ttɛawanen asekyed nniḍen; amussnaw isbadu akayad nniḍen i ilaqen.

Iɣbula d yisemliyen

Ilaq uqcic ad yeqqim iman-is?

Ameskar isbadu aktil seg wazal n 1 m, ula i lufan i yeṭṭef win it-iɛawnen. Amussnaw iheggi asers d uṭṭaf n twali i unṭaf.

Iɣbula d yisemliyen

Zemren yigmaḍ ad ttusekyeden s PC?

Asebter n umeskar isbadu asekyed s PC akked HIIS-1 d umasal DICOM. Aseɣẓan d uẓeṭṭa ilaq ad ttuseggmen i umecwar-a.

Iɣbula d yisemliyen

Wali akatalug akliniki n ASEM

Wali taɣulin icudden deg ukatalug ASEM. Iseɣwan-a ur rennun ara isnasen imaynuten i yibenk-a.

Asenqed d biyumetri

Ugar n yibenk, d amɛiwen i wakud ɣezzifen

Mmeslay ɣef usenfar-ik akked ASEM: afran n yibenkan, asbeddi, asileɣ d tallalt tatiknikant.



Tugna n umeẓlu

Asbeddi

Mmeslay ɣef yisuturen n usbeddi d twila akked terbaɛt-nneɣ tatiknikant.

Issin ugar



Tugna n umeẓlu

Asileɣ

Suter isallen ɣef usileɣ d ulmad n useqdec yellan i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Ameẓlu deffir uznuzu

Nermes ASEM i tallalt tatiknikant d umɛiwen i yibenk-ik.

Issin ugar



Tugna n umeẓlu

Anṣaḥ n yimussnawen

Mmeslay-aɣ ɣef yisuturen-ik iklinikiyen i wakken ad nefren ibenkan iwulmen.

Issin ugar

Nermes amussnaw

Wali akk ifarisen

Asmel n umeskar

+21320092281 · direction@asem-medical.com

Lqem n yisemli: 13414acc5c44